Package 'ProData'

December 4, 2025

,	
Title SELDI-TOF data of Breast cancer samples	
Version 1.48.0	
Author Xiaochun Li	
Description A data package of SELDI-TOF protein mass spectrometry data of 167 breast cancer and normal samples.	
Maintainer Xiaochun Li <xiaochun@jimmy.harvard.edu></xiaochun@jimmy.harvard.edu>	
Depends R (>= 2.4.0), Biobase (>= 2.5.5)	
License GPL	
biocViews ExperimentData, CancerData, BreastCancerData, MassSpectrometryData, NCI	
git_url https://git.bioconductor.org/packages/ProData	
git_branch RELEASE_3_22	
git_last_commit 141570f git_last_commit_date 2025-10-29 Repository Bioconductor 3.22 Date/Publication 2025-12-04	
Contents	
f45cbmk	
Index	4
f 45cbmk A SELDI-TOF data of 167 breast cancer and normal samples from a Miron Laboratory	the

2 ProData

Description

The data consist of 167 SELDI-TOF spectra from either individuals with breast cancer or without (normal controls). Two covariates are available. The data have been pre-processed by Ciphergen's software and the specific details are as follows. The baselines of the spectra were subtracted using the default parameter settings of the Ciphergen software, and the spectra were normalized by the total ion current normalization function of the Ciphergen software. Peaks were detected with signal to noise ratio (S/N) above 10. Proto-biomarkers, i.e., peaks in each spectrum that have the same M/Z value, were detected with the following parameters: first pass, 10 S/N, min peak threshold 10 pass, 5 S/N, and add estimated peaks to complete clusters and it is the jointly normalized data that are available here. The data are presented in the form of an exprSet object.

Usage

data(f45cbmk)

Format

The different covariates are:

GROUP a character vector with four possible values, "A", "B", "C" and "D". "A" represents HER2 positive patients, "B" normal healthy patients, "C" ER/PR positive patients and "D" samples from a single healthy woman.

SPECTRUM a character vector, the spectrum id of each sample.

Source

Alex Miron's Lab at the Dana-Farber Cancer Institute.

References

Shi, Q. et al. Biomarker Discovery in Plasma Using Proteomics Analysis for Early Detection of Breast Cancer. 2004, manuscript.

Examples

data(f45cbmk)

ProData

A SELDI-TOF data of 167 breast cancer and normal spectra

Description

This library contains both the raw spectra and the proto-biomarker data pre-processed by Ciphergen's software of the 167 samples. Those one hundred and sixty seven samples are collected from 155 subjects in CPT tubes with plasma isolated and stored in -80C until needed. Among the 167 samples, 55 are HER2 positive (A), 64 are normal healthy women (B), 35 are ER/PR positive (mostly) (C) and 13 samples are from a single healthy woman. Samples in group D are the only ones from a single subject, all the other samples represent draws from individual subjects. Samples were thawed and aliquoted into 100ul vials. The samples were fractionated to simplify the proteome into sub-proteomes. The fractions 4 and 5 (f45) were processed by the Ciphergen IMAC protocol with EAM of CHCA.

ProData 3

Information on the spectrum ID and the pheotype information is stored in the exprSet object "f45cbmk", the proto-biomakers pre-processed by Ciphergen's software.

As an alternative, package PROcess may be used for pre-processing of the raw spectra to get a set of proto-biomarkers.

Source

Alex Miron's Lab at the Dana-Farber Cancer Institute.

Examples

```
# plot a raw spectrum
f45c <- system.file("f45c", package="ProData")
fs <- dir(f45c,full.names=TRUE)
plot(read.csv(gzfile(fs[1])), type="1")

# find out sizes of phenotype groups:
library(Biobase)
data(f45cbmk)
SpecGrp <- pData(f45cbmk)
table(SpecGrp[,1])</pre>
```

Index

```
* datasets
f45cbmk, 1
ProData, 2
f45cbmk, 1
ProData, 2
```