

# Package ‘CytoPipeline’

December 4, 2025

**Title** Automation and visualization of flow cytometry data analysis pipelines

**Version** 1.11.0

**Description** This package provides support for automation and visualization of flow cytometry data analysis pipelines. In the current state, the package focuses on the preprocessing and quality control part. The framework is based on two main S4 classes, i.e. CytoPipeline and CytoProcessingStep. The pipeline steps are linked to corresponding R functions - that are either provided in the CytoPipeline package itself, or exported from a third party package, or coded by the user her/himself. The processing steps need to be specified centrally and explicitly using either a json input file or through step by step creation of a CytoPipeline object with dedicated methods. After having run the pipeline, obtained results at all steps can be retrieved and visualized thanks to file caching (the running facility uses a BiocFileCache implementation). The package provides also specific visualization tools like pipeline workflow summary display, and 1D/2D comparison plots of obtained flowFrames at various steps of the pipeline.

**License** GPL-3

**Encoding** UTF-8

**Roxygen** list(markdown = TRUE)

**RoxygenNote** 7.3.2

**BugReports** <https://github.com/UCLouvain-CBIO/CytoPipeline/issues>

**URL** <https://uclouvain-cbio.github.io/CytoPipeline>

**biocViews** FlowCytometry, Preprocessing, QualityControl, WorkflowStep, ImmunoOncology, Software, Visualization

**Collate** 'CytoPipeline-functions.R' 'CytoPipeline-package.R'  
'CytoPipelineClass.R' 'CytoProcessingStep.R'  
'CytoProcessingStepImplementations.R' 'data.R' 'gating.R'  
'utils.R' 'ggplots.R'

**Depends** R (>= 4.4)

**Imports** methods, stats, utils, withr, rlang, ggplot2 (>= 3.4.1),  
ggcyto, BiocFileCache, BiocParallel, flowCore, PeacoQC, flowAI,  
diagram, jsonlite, scales

**Suggests** testthat (>= 3.0.0), vdiff, diffviewer, knitr, rmarkdown,  
BiocStyle, reshape2, dplyr, CytoPipelineGUI

**VignetteBuilder** knitr

**Config/testthat/edition** 3

**git\_url** <https://git.bioconductor.org/packages/CytoPipeline>

**git\_branch** devel

**git\_last\_commit** 54ed268

**git\_last\_commit\_date** 2025-10-29

**Repository** Bioconductor 3.23

**Date/Publication** 2025-12-04

**Author** Philippe Hauchamps [aut, cre] (ORCID:  
<<https://orcid.org/0000-0003-2865-1852>>),  
Laurent Gatto [aut] (ORCID: <<https://orcid.org/0000-0002-1520-2268>>),  
Dan Lin [ctb]

**Maintainer** Philippe Hauchamps <[philippe.hauchamps@uclouvain.be](mailto:philippe.hauchamps@uclouvain.be)>

## Contents

|                                             |    |
|---------------------------------------------|----|
| aggregateAndSample . . . . .                | 3  |
| appendCellID . . . . .                      | 4  |
| applyScaleTransforms . . . . .              | 5  |
| areFluoCols . . . . .                       | 6  |
| areSignalCols . . . . .                     | 6  |
| compensateFromMatrix . . . . .              | 7  |
| computeScatterChannelsLinearScale . . . . . | 8  |
| CytoPipeline . . . . .                      | 9  |
| CytoPipeline-class . . . . .                | 10 |
| CytoProcessingStep . . . . .                | 16 |
| estimateScaleTransforms . . . . .           | 17 |
| execute . . . . .                           | 19 |
| exportCytoPipeline . . . . .                | 24 |
| findTimeChannel . . . . .                   | 25 |
| getAcquiredCompensationMatrix . . . . .     | 26 |
| getChannelNamesFromMarkers . . . . .        | 27 |
| getFCSFileName . . . . .                    | 28 |
| getTransfoParams . . . . .                  | 28 |
| ggplotEvents . . . . .                      | 30 |
| ggplotFilterEvents . . . . .                | 34 |
| ggplotFlowRate . . . . .                    | 36 |
| handlingProcessingSteps . . . . .           | 37 |
| inspectCytoPipelineObjects . . . . .        | 40 |
| interactingWithCytoPipelineCache . . . . .  | 45 |
| OMIP021Samples . . . . .                    | 47 |
| qualityControlFlowAI . . . . .              | 47 |
| qualityControlPeacoQC . . . . .             | 49 |

|                                      |    |
|--------------------------------------|----|
| readRDSObject . . . . .              | 50 |
| readSampleFiles . . . . .            | 51 |
| removeChannels . . . . .             | 53 |
| removeDeadCellsManualGate . . . . .  | 53 |
| removeDebrisManualGate . . . . .     | 55 |
| removeDoubletsCytoPipeline . . . . . | 56 |
| removeMarginsPeacoQC . . . . .       | 57 |
| resetCellIDs . . . . .               | 58 |
| runCompensation . . . . .            | 59 |
| singletsGate . . . . .               | 60 |
| subsample . . . . .                  | 62 |
| updateMarkerName . . . . .           | 63 |
| writeFlowFrame . . . . .             | 63 |

|              |           |
|--------------|-----------|
| <b>Index</b> | <b>66</b> |
|--------------|-----------|

---

|                    |                                                                         |
|--------------------|-------------------------------------------------------------------------|
| aggregateAndSample | <i>Aggregate and sample multiple flow frames of a flow set together</i> |
|--------------------|-------------------------------------------------------------------------|

---

## Description

Aggregate multiple flow frames in order to analyze them simultaneously. A new FF, which contains about  $nTotalEvents$  cells,  $nTotalEvents/nFiles$  cells from each file. Two new columns are added: a column indicating the original file by index, and a noisy version of this, for better plotting opportunities. This function is based on `PeacoQC::AggregateFlowframes()` where file names inputs have been replaced by a `flowSet` input.

## Usage

```
aggregateAndSample(
  fs,
  nTotalEvents,
  setup = c("forceNEvent", "forceBalance"),
  seed = NULL,
  channels = NULL,
  writeOutput = FALSE,
  outputFile = "aggregate.fcs",
  keepOrder = FALSE
)
```

## Arguments

|                           |                                                                                    |
|---------------------------|------------------------------------------------------------------------------------|
| <code>fs</code>           | a <code>flowCore::flowset</code>                                                   |
| <code>nTotalEvents</code> | Total number of cells to select from the input flow frames                         |
| <code>setup</code>        | How to proceed when $nTotalEvents/nFiles$ is too high for some of the flow frames: |

|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                    | <ul style="list-style-type: none"> <li>• <b>forceBalance</b> (default): compute the minimum nb of events per flow frame, and keep that amount of events from each flow frame.</li> <li>• <b>forceNEvents</b>: try to be as balanced as possible, but force a total of <b>nTotalEvents</b> if possible, i.e. takes all events from the flow frame with too low nb of events, and then fill in the total with events from the bigger flow frames in a balanced way. However, if <b>nTotalEvents</b> is greater than the sum of all events, take all events only once.</li> </ul> |
| <b>seed</b>        | seed to be set before sampling for reproducibility. Default NULL does not set any seed.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>channels</b>    | Channels/markers to keep in the aggregate. Default NULL takes all channels of the first file.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>writeOutput</b> | Whether to write the resulting flowframe to a file. Default FALSE                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>outputFile</b>  | Full path to output file. Default "aggregate.fcs"                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>keepOrder</b>   | If TRUE, the random subsample will be ordered in the same way as they were originally ordered in the file. Default = FALSE.                                                                                                                                                                                                                                                                                                                                                                                                                                                    |

### Value

returns a new `flowCore::flowFrame`

### Examples

```
data(OMIP021Samples)

nCells <- 1000
agg <- aggregateAndSample(
  fs = OMIP021Samples,
  nTotalEvents = nCells)
```

---

|                           |                                                   |
|---------------------------|---------------------------------------------------|
| <code>appendCellID</code> | <i>append 'Original_ID' column to a flowframe</i> |
|---------------------------|---------------------------------------------------|

---

### Description

: on a `flowCore::flowFrame`, append a 'Original\_ID' column. This column can be used in plots comparing the events pre and post gating. If the 'Original\_ID' column already exists, the function does nothing

### Usage

```
appendCellID(ff, eventIDs = seq_len(flowCore::nrow(ff)))
```

### Arguments

|                 |                                                                                             |
|-----------------|---------------------------------------------------------------------------------------------|
| <b>ff</b>       | a <code>flowCore::flowFrame</code>                                                          |
| <b>eventIDs</b> | an integer vector containing the values to be added in expression matrix, as Original ID's. |

## Value

new flowCore::flowFrame containing the added 'Original\_ID' column

## Examples

```
data(OMIP021Samples)

retFF <- appendCellID(OMIP021Samples[[1]])
```

---

applyScaleTransforms    *apply scale transforms*

---

## Description

wrapper around flowCore::transform() that discards any additional parameter passed in (...) Additionally, some checks regarding channels correspondance is done: if transList contains transformations for channels that are not present in x, then these transformations are first removed.

## Usage

```
applyScaleTransforms(x, transList, verbose = FALSE, ...)
```

## Arguments

|           |                                                   |
|-----------|---------------------------------------------------|
| x         | a flowCore::flowSet or a flowCore::flowFrame      |
| transList | a flowCore::transformList                         |
| verbose   | if TRUE, send a message per flowFrame transformed |
| ...       | other arguments (not used)                        |

## Value

the transformed flowFrame

## Examples

```
data(OMIP021Samples)

transListPath <- file.path(system.file("extdata",
                                         package = "CytoPipeline"),
                           "OMIP021_TransList.rds")

transList <- readRDSObject(transListPath)

ff_c <- compensateFromMatrix(OMIP021Samples[[1]],
                             matrixSource = "fcs")

ff_t <- applyScaleTransforms(ff_c, transList = transList)
```

---

|             |                                                                    |
|-------------|--------------------------------------------------------------------|
| areFluoCols | <i>find flow frame columns that represent fluorochrome channel</i> |
|-------------|--------------------------------------------------------------------|

---

**Description**

: find flow frame columns that represent fluorochrome channel

**Usage**

```
areFluoCols(
  x,
  toRemovePatterns = c("FSC", "SSC", "Time", "Original_ID", "File", "SampleID")
)
```

**Arguments**

**x** a flowCore::flowFrame or a flowCore::flowSet  
**toRemovePatterns** a vector of string patterns that are to be considered as non fluorochrome

**Value**

a vector of booleans of which the dimension is equal to the number of columns in ff

**Examples**

```
data(OMIP021Samples)

areFluoCols(OMIP021Samples)
```

---

|               |                                                           |
|---------------|-----------------------------------------------------------|
| areSignalCols | <i>find flow frame columns that represent true signal</i> |
|---------------|-----------------------------------------------------------|

---

**Description**

: find flow frame columns that represent true signal

**Usage**

```
areSignalCols(
  x,
  toRemovePatterns = c("Time", "Original_ID", "File", "SampleID")
)
```

**Arguments**

`x` a `flowCore::flowFrame` or a `flowCore::flowSet`  
`toRemovePatterns` a vector of string patterns that are to be considered as non signal

**Value**

a vector of booleans of which the dimension is equal to the number of columns in `ff`

**Examples**

```
data(OMIP021Samples)
areSignalCols(OMIP021Samples)
```

---

`compensateFromMatrix`    *compensation of fcs file(s) from matrix*

---

**Description**

executes the classical compensation function on a `flowSet` or `flowFrame`, given a compensation matrix. The matrix can be either retrieved in the fcs files themselves or provided as a csv file.

**Usage**

```
compensateFromMatrix(
  x,
  matrixSource = c("fcs", "import"),
  matrixPath = NULL,
  updateChannelNames = TRUE,
  verbose = FALSE,
  ...
)
```

**Arguments**

`x` a `flowCore::flowFrame` or `flowCore::flowSet`  
`matrixSource` if "fcs", the compensation matrix will be fetched from the fcs files (different compensation matrices can then be applied by fcs file) if "import", uses `matrixPath` to read the matrix (should be a csv file)  
`matrixPath` if `matrixSource == "import"`, will be used as the input csv file path  
`updateChannelNames` if TRUE, updates the fluo channel names by prefixing them with "comp-"  
`verbose` if TRUE, displays information messages  
`...` additional arguments (not used)

**Value**

the compensated flowSet or flowFrame

**Examples**

```
rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
sampleFiles <-
  file.path(rawDataDir, list.files(rawDataDir, pattern = "Donor"))

truncateMaxRange <- FALSE
minLimit <- NULL

# create flowCore::flowSet with all samples of a dataset
fsRaw <- readSampleFiles(
  sampleFiles = sampleFiles,
  whichSamples = "all",
  truncate_max_range = truncateMaxRange,
  min.limit = minLimit)

suppressWarnings(ff_m <- removeMarginsPeacoQC(x = fsRaw[[2]]))

ff_c <-
  compensateFromMatrix(ff_m,
    matrixSource = "fcs")
```

---

**computeScatterChannelsLinearScale**

*compute linear transformation of scatter channels found in ff, based on 5% and 95% of referenceChannel, set as target. If there is a transformation defined in transList for referenceChannel, it is applied first, before computing quantiles. Then the computed linear transformations (or each scatter channel) are added into the transfo\_list. -A channels are computed, and same linear transformation is then applied to corresponding -W and -H channels (if they exist in ff).*

---

**Description**

based on a referenceChannel

**Usage**

```
computeScatterChannelsLinearScale(
  ff,
  transList = NULL,
  referenceChannel,
  silent = TRUE
)
```

**Arguments**

|                               |                                                                                                      |
|-------------------------------|------------------------------------------------------------------------------------------------------|
| <code>ff</code>               | a <code>flowCore::flowFrame</code>                                                                   |
| <code>transList</code>        | an initial <code>flowCore::transformList</code>                                                      |
| <code>referenceChannel</code> | the reference channel to take target quantile values from. Can be defined as marker or channel name. |
| <code>silent</code>           | if FALSE, will output some information on the computed linear transformations                        |

**Value**

the `transList` with added linear scale transformations

**Examples**

```
data(OMIP021Samples)

ff <- OMIP021Samples[[1]]
refMarker <- "APCCy7 - CD4"
refChannel <- "780/60Red-A"
transList <- flowCore::estimateLogicle(ff,
                                     channels = refChannel)

retTransList <-
  computeScatterChannelsLinearScale(ff,
                                   transList = transList,
                                   referenceChannel = refMarker,
                                   silent = TRUE
  )
```

---

CytoPipeline

*CytoPipeline package*

---

**Description**

CytoPipeline is a package that provides support for automation and visualization of flow cytometry data analysis pipelines. In the current state, the package focuses on the preprocessing and quality control part.

The framework is based on two main S4 classes, i.e. `CytoPipeline` and `CytoProcessingStep`. The `CytoProcessingStep` defines the link between pipeline step names and corresponding R functions that are either provided in the `CytoPipeline` package itself, or exported from a third party package, or coded by the user her/himself. The processing steps need to be specified centrally and explicitly using either a json input file or through step by step creation of a `CytoPipeline` object with dedicated methods.

After having run the pipeline, obtained results at all steps can be retrieved and visualized thanks to file caching (the running facility uses a `BiocFileCache` implementation). The package provides also specific visualization tools like pipeline workflow summary display, and 1D/2D comparison plots of obtained `flowFrames` at various steps of the pipeline.

For a step by step example using `CytoPipeline`, please have a look at the vignette!

**Author(s)**

**Maintainer:** Philippe Hauchamps <philippe.hauchamps@uclouvain.be> ([ORCID](#))

Authors:

- Laurent Gatto <laurent.gatto@uclouvain.be> ([ORCID](#))

Other contributors:

- Dan Lin <dan.8.lin@gsk.com> [contributor]

**See Also**

[CytoPipelineClass](#), [CytoProcessingStep](#)

---

|                    |                           |
|--------------------|---------------------------|
| CytoPipeline-class | <i>CytoPipeline class</i> |
|--------------------|---------------------------|

---

**Description**

Class representing a flow cytometry pipeline, and composed of two processing queues, i.e. lists of CytoProcessingStep objects :

- a list of CytoProcessingStep(s) for pre-calculation of scale transformations per channel
- a list of CytoProcessingStep(s) for the pre-processing of flow frames

**Usage**

```
## S4 method for signature 'CytoPipeline'
show(object)

## S4 method for signature 'missing'
CytoPipeline(
  object,
  experimentName = "default_experiment",
  sampleFiles = character(),
  pData = NULL
)

## S4 method for signature 'list'
CytoPipeline(
  object,
  experimentName = "default_experiment",
  sampleFiles = character(),
  pData = NULL
)

## S4 method for signature 'character'
```

```

CytoPipeline(
  object,
  experimentName = "default_experiment",
  sampleFiles = character(),
  pData = NULL
)

## S3 method for class 'CytoPipeline'
as.list(x, ...)

experimentName(x)

experimentName(x) <- value

sampleFiles(x)

sampleFiles(x) <- value

pData(x)

pData(x) <- value

sampleDisplayNames(x, sampleFiles = NULL)

sampleNameFromDisplayName(x, displayName)

```

### Arguments

|                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| object         | a character() containing a JSON input                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| experimentName | the experiment name                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| sampleFiles    | a character (e.g. sampleFileNames) or a numeric vector (e.g. indices of sample files). If NULL, all samples will be displayed.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| pData          | the pheno Data (data.frame or NULL)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| x              | a CytoPipeline object                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| ...            | additional arguments (not used here)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| value          | the new value to be assigned. the pData<- setter is a bit more liberal than it used to be: <ol style="list-style-type: none"> <li>1. It can accept new pData containing more rows than existing sample names (the corresponding subset of pData is taken).</li> <li>2. It can accept pData with row names pointing to either sample file full paths or base file names</li> <li>3. It can accept pData with no row names provided the number of rows correspond to the number of sample files. Row names are then set by default to sample file base names (if unique), or sample file full paths.</li> </ol> |
| displayName    | a character                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |

**Value**

nothing

- for `as.list.CytoPipeline`: the obtained list
- for `sampleDisplayNames`: a character vector of sample display names
- for `sampleNameFromDisplayName`: the sample name corresponding to the specified display name. of sample display names

**Slots**

`scaleTransformProcessingQueue` A list of `CytoProcessingStep` objects containing the steps for obtaining the scale transformations per channel

`flowFramesPreProcessingQueue` A list of `CytoProcessingStep` objects containing the steps for pre-processing of the samples flow frames

`experimentName` A character containing the experiment (run) name

`sampleFiles` A character vector storing all fcs files to be run into the pipeline

`pData` An optional data.frame containing additional information for each sample file. The `pData` raw names should correspond to the sample files (using full paths or base paths). If the `pData` contains a columns with name 'displayName', this will have an impact in the `sampleDisplayNames()` function, i.e. sample display names will be the one mentioned in `pData`, instead of typically base file names (or larger paths if base file names are not unique)

**Examples**

```
### *** EXAMPLE 1: building CytoPipeline step by step *** ###

rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
experimentName <- "OMIP021_PeacoQC"
sampleFiles <- file.path(rawDataDir, list.files(rawDataDir,
                                                pattern = "Donor"))

outputDir <- base::tempdir()

# main parameters : sample files and output files
pipL <- CytoPipeline(experimentName = experimentName,
                    sampleFiles = sampleFiles)

### SCALE TRANSFORMATION STEPS ###

pipL <-
  addProcessingStep(pipL,
                    whichQueue = "scale transform",
                    CytoProcessingStep(
                      name = "flowframe_read",
                      FUN = "readSampleFiles",
                      ARGS = list(
                        whichSamples = "all",
```

```

        truncate_max_range = FALSE,
        min.limit = NULL
      )
    )
  )

  pipl <-
    addProcessingStep(pipl,
      whichQueue = "scale transform",
      CytoProcessingStep(
        name = "remove_margins",
        FUN = "removeMarginsPeacoQC",
        ARGS = list()
      )
    )

  pipl <-
    addProcessingStep(pipl,
      whichQueue = "scale transform",
      CytoProcessingStep(
        name = "compensate",
        FUN = "compensateFromMatrix",
        ARGS = list(matrixSource = "fcs")
      )
    )

  pipl <-
    addProcessingStep(pipl,
      whichQueue = "scale transform",
      CytoProcessingStep(
        name = "flowframe_aggregate",
        FUN = "aggregateAndSample",
        ARGS = list(
          nTotalEvents = 10000,
          seed = 0
        )
      )
    )

  pipl <-
    addProcessingStep(pipl,
      whichQueue = "scale transform",
      CytoProcessingStep(
        name = "scale_transform_estimate",
        FUN = "estimateScaleTransforms",
        ARGS = list(
          fluoMethod = "estimateLogicle",
          scatterMethod = "linear",
          scatterRefMarker = "BV785 - CD3"
        )
      )
    )

```

```

### PRE-PROCESSING STEPS ###

piPL <-
  addProcessingStep(piPL,
    whichQueue = "pre-processing",
    CytoProcessingStep(
      name = "flowframe_read",
      FUN = "readSampleFiles",
      ARGS = list(
        truncate_max_range = FALSE,
        min.limit = NULL
      )
    )
  )

piPL <-
  addProcessingStep(piPL,
    whichQueue = "pre-processing",
    CytoProcessingStep(
      name = "remove_margins",
      FUN = "removeMarginsPeacoQC",
      ARGS = list()
    )
  )

piPL <-
  addProcessingStep(piPL,
    whichQueue = "pre-processing",
    CytoProcessingStep(
      name = "compensate",
      FUN = "compensateFromMatrix",
      ARGS = list(matrixSource = "fcs")
    )
  )

piPL <-
  addProcessingStep(
    piPL,
    whichQueue = "pre-processing",
    CytoProcessingStep(
      name = "remove_debris",
      FUN = "removeDebrisManualGate",
      ARGS = list(
        FSCChannel = "FSC-A",
        SSCChannel = "SSC-A",
        gateData = c(73615, 110174, 213000, 201000, 126000,
                     47679, 260500, 260500, 113000, 35000)))
  )

piPL <-
  addProcessingStep(piPL,
    whichQueue = "pre-processing",
    CytoProcessingStep(

```

```

        name = "remove_dead_cells",
        FUN = "removeDeadCellsManualGate",
        ARGS = list(
            FSCChannel = "FSC-A",
            LDMarker = "L/D Aqua - Viability",
            gateData = c(0, 0, 250000, 250000,
                        0, 650, 650, 0)
        )
    )
)

pipl <-
  addProcessingStep(
    pipl,
    whichQueue = "pre-processing",
    CytoProcessingStep(
      name = "perform_QC",
      FUN = "qualityControlPeacoQC",
      ARGS = list(
        preTransform = TRUE,
        min_cells = 150, # default
        max_bins = 500, # default
        step = 500, # default,
        MAD = 6, # default
        IT_limit = 0.55, # default
        force_IT = 150, # default
        peak_removal = 0.3333, # default
        min_nr_bins_peakdetection = 10 # default
      )
    )
  )

pipl <-
  addProcessingStep(pipl,
    whichQueue = "pre-processing",
    CytoProcessingStep(
      name = "transform",
      FUN = "applyScaleTransforms",
      ARGS = list()
    )
  )

### *** EXAMPLE 2: building CytoPipeline from JSON file *** ###

jsonDir <- system.file("extdata", package = "CytoPipeline")
jsonPath <- file.path(jsonDir, "pipelineParams.json")

pipL2 <- CytoPipeline(jsonPath,
  experimentName = experimentName,
  sampleFiles = sampleFiles)

```

---

|                    |                             |
|--------------------|-----------------------------|
| CytoProcessingStep | <i>Cyto Processing step</i> |
|--------------------|-----------------------------|

---

## Description

Class containing the function and arguments to be applied in a lazy-execution framework.

Objects of this class are created using the `CytoProcessingStep()` function. The processing step is executed with the `executeProcessingStep()` function.

## Usage

```
CytoProcessingStep(name = character(), FUN = character(), ARGS = list())
```

```
## S4 method for signature 'CytoProcessingStep'
show(object)
```

```
executeProcessingStep(x, ...)
```

```
getCPSName(x)
```

```
getCPSFUN(x)
```

```
getCPSARGS(x)
```

```
## S3 method for class 'CytoProcessingStep'
as.list(x, ...)
```

```
as.json.CytoProcessingStep(x, pretty = FALSE)
```

```
from.json.CytoProcessingStep(jsonString)
```

## Arguments

|            |                                                                                      |
|------------|--------------------------------------------------------------------------------------|
| name       | character denoting a name to the step, which can be different from the function name |
| FUN        | function or character representing a function name.                                  |
| ARGS       | list of arguments to be passed along to FUN.                                         |
| object     | a <code>CytoProcessingStep</code> object.                                            |
| x          | a <code>CytoProcessingStep</code> object.                                            |
| ...        | other arguments (not used)                                                           |
| pretty     | formatting set-up (see <code>jsonlite::toJSON</code> doc)                            |
| jsonString | a <code>character()</code> containing a JSON string.                                 |

## Details

This object contains all relevant information of a data analysis processing step, i.e. the function and all of its arguments to be applied to the data.

## Value

The CytoProcessingStep function returns an object of type CytoProcessingStep.

## Examples

```
## Create a simple processing step object
ps1 <- CytoProcessingStep("summing step", sum)

getCPSName(ps1)

getCPSFUN(ps1)

getCPSARGS(ps1)

executeProcessingStep(ps1, 1:10)

as.list(ps1)

js_str <- as.json.CytoProcessingStep(ps1)

ps2 <- from.json.CytoProcessingStep(js_str)

identical(ps1, ps2)
```

---

estimateScaleTransforms

*estimates scale transformations*

---

## Description

this function estimates the scale transformations to be applied on a flowFrame to obtain 'good behaving' distributions, i.e. the best possible separation between + population and - population. It distinguishes between scatter channels, where either linear, or no transform is applied, and fluo channels, where either logicle transform

- using flowCore::estimateLogicle - is estimated, or no transform is applied.

The idea of linear transform of scatter channels is as follows: a reference channel (not a scatter one) is selected and a linear transform ( $Y = AX + B$ ) is applied to all scatter channel, as to align their 5 and 95 percentiles to those of the reference channel. For the estimateLogicle function, see flowCore documentation.

**Usage**

```
estimateScaleTransforms(
  ff,
  fluoMethod = c("estimateLogicle", "none"),
  scatterMethod = c("none", "linearQuantile"),
  scatterRefMarker = NULL,
  specificScatterChannels = NULL,
  verbose = FALSE
)
```

**Arguments**

|                                      |                                                                                                          |
|--------------------------------------|----------------------------------------------------------------------------------------------------------|
| <code>ff</code>                      | a <code>flowCore::flowFrame</code>                                                                       |
| <code>fluoMethod</code>              | method to be applied to all fluo channels                                                                |
| <code>scatterMethod</code>           | method to be applied to all scatter channels                                                             |
| <code>scatterRefMarker</code>        | the reference channel that is used to align the                                                          |
| <code>specificScatterChannels</code> | vector of scatter channels for which we still want to apply the fluo method (and not the scatter Method) |
| <code>verbose</code>                 | if TRUE, send messages to the user at each step                                                          |

**Value**

a `flowCore::flowFrame` with removed low quality events from the input

**Examples**

```
data(OMIP021Samples)

compMatrix <- flowCore::spillover(OMIP021Samples[[1]])$SPILL
ff_c <- runCompensation(OMIP021Samples[[1]], spillover = compMatrix)

transList <-
  estimateScaleTransforms(
    ff = ff_c,
    fluoMethod = "estimateLogicle",
    scatterMethod = "linear",
    scatterRefMarker = "BV785 - CD3")
```

---

|         |                                      |
|---------|--------------------------------------|
| execute | <i>executing CytoPipeline object</i> |
|---------|--------------------------------------|

---

## Description

this function triggers the execution of the processing queues of a CytoPipeline object. First, the scale tranform processing queue is run, taking the set of sample names as an implicit first input. At the end of the queue, a scale transform List is assumed to be created. Second, the flowFrame pre-processing queue, repeatedly for each sample file. The scale transform list generated in the previous step is taken as implicit input, together with the initial sample file. At the end of the queue run, a pre-processed flowFrame is assumed to be generated. No change is made on the input CytoPipeline object, all results are stored in the cache.

## Usage

```
execute(
  x,
  path = ".",
  rmCache = FALSE,
  useBiocParallel = FALSE,
  BPPARAM = BiocParallel::bpparam(),
  BPOPTIONS = BiocParallel::bpoptions(packages = c("flowCore")),
  saveLastStepFF = TRUE,
  saveFFSuffix = "_preprocessed",
  saveFFFormat = c("fcs", "csv"),
  saveFFCsvUseChannelMarker = TRUE,
  saveScaleTransforms = FALSE
)
```

## Arguments

|                 |                                                                                                                                                                            |
|-----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| x               | CytoPipeline object                                                                                                                                                        |
| path            | base path, a subdirectory with name equal to the experiment will be created to store the output data, in particular the experiment cache                                   |
| rmCache         | if TRUE, starts by removing the already existing cache directory corresponding to the experiment                                                                           |
| useBiocParallel | if TRUE, use BiocParallel for computation of the sample file pre-processing in parallel (one file per worker at a time). Note the BiocParallel function used is bplapply() |
| BPPARAM         | if useBiocParallel is TRUE, sets the BPPARAM back-end to be used for the computation. If not provided, will use the top back-end on the BiocParallel::registered() stack.  |
| BPOPTIONS       | if useBiocParallel is TRUE, sets the BPOPTIONS to be passed to bplapply() function. Note that if you use a SnowParams back-end, you need to specify                        |

all the packages that need to be loaded for the different `CytoProcessingStep` to work properly (visibility of functions). As a minimum, the `flowCore` package needs to be loaded. (hence the default `BPOPTIONS = bpoptions(packages = c("flowCore"))`)

`saveLastStepFF` if TRUE, save the final result of the pre-processing, for each file. By convention, these output files are stored in `path/x@experimentName/output/`, the file names used are the same as the initial fcs file basenames, concatenated with `saveFFSuffix`, and with file extension corresponding to `saveFFFormat`.

`saveFFSuffix` FF file name suffix

`saveFFFormat` either fcs or csv

`saveFFCsvUseChannelMarker` if TRUE (default), converts the channels to the corresponding marker names (where the Marker is not NA). This setting is only applicable to export in csv format.

`saveScaleTransforms` if TRUE (default FALSE), save on disk (in RDS format) the `flowCore::transformList` object obtained after running the `scaleTransform` processing queue. The file name is hardcoded to `path/experimentName/RDS/scaleTransformList.rds`

## Value

nothing

## Examples

```
### *** EXAMPLE 1: building CytoPipeline step by step *** ###

rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
experimentName <- "OMIP021_PeacoQC"
sampleFiles <- file.path(rawDataDir, list.files(rawDataDir,
  pattern = "Donor"))

outputDir <- base::tempdir()

# main parameters : sample files and output files
pipelineParams <- list()
pipelineParams$experimentName <- experimentName
pipelineParams$sampleFiles <- sampleFiles
pipl <- CytoPipeline(pipelineParams)

### SCALE TRANSFORMATION STEPS ###

pipl <-
  addProcessingStep(pipl,
    whichQueue = "scale transform",
    CytoProcessingStep(
      name = "flowframe_read",
      FUN = "readSampleFiles",
      ARGS = list(
```

```

        whichSamples = "all",
        truncate_max_range = FALSE,
        min.limit = NULL
    )
)

)

pipL <-
  addProcessingStep(pipL,
    whichQueue = "scale transform",
    CytoProcessingStep(
      name = "remove_margins",
      FUN = "removeMarginsPeacoQC",
      ARGS = list()
    )
  )

pipL <-
  addProcessingStep(pipL,
    whichQueue = "scale transform",
    CytoProcessingStep(
      name = "compensate",
      FUN = "compensateFromMatrix",
      ARGS = list(matrixSource = "fcs")
    )
  )

pipL <-
  addProcessingStep(pipL,
    whichQueue = "scale transform",
    CytoProcessingStep(
      name = "flowframe_aggregate",
      FUN = "aggregateAndSample",
      ARGS = list(
        nTotalEvents = 10000,
        seed = 0
      )
    )
  )

pipL <-
  addProcessingStep(pipL,
    whichQueue = "scale transform",
    CytoProcessingStep(
      name = "scale_transform_estimate",
      FUN = "estimateScaleTransforms",
      ARGS = list(
        fluoMethod = "estimateLogicle",
        scatterMethod = "linear",
        scatterRefMarker = "BV785 - CD3"
      )
    )
  )

```

```

### PRE-PROCESSING STEPS ###

pipl <-
  addProcessingStep(pipl,
    whichQueue = "pre-processing",
    CytoProcessingStep(
      name = "flowframe_read",
      FUN = "readSampleFiles",
      ARGS = list(
        truncate_max_range = FALSE,
        min.limit = NULL
      )
    )
  )

pipl <-
  addProcessingStep(pipl,
    whichQueue = "pre-processing",
    CytoProcessingStep(
      name = "remove_margins",
      FUN = "removeMarginsPeacoQC",
      ARGS = list()
    )
  )

pipl <-
  addProcessingStep(pipl,
    whichQueue = "pre-processing",
    CytoProcessingStep(
      name = "compensate",
      FUN = "compensateFromMatrix",
      ARGS = list(matrixSource = "fcs")
    )
  )

pipl <-
  addProcessingStep(
    pipl,
    whichQueue = "pre-processing",
    CytoProcessingStep(
      name = "remove_debris",
      FUN = "removeDebrisManualGate",
      ARGS = list(
        FSCChannel = "FSC-A",
        SSCChannel = "SSC-A",
        gateData = c(73615, 110174, 213000, 201000, 126000,
                     47679, 260500, 260500, 113000, 35000)
      )
    )
  )

pipl <-

```

```

        addProcessingStep(pipL,
                          whichQueue = "pre-processing",
                          CytoProcessingStep(
                            name = "remove_dead_cells",
                            FUN = "removeDeadCellsManualGate",
                            ARGS = list(
                              FSCChannel = "FSC-A",
                              LDMarker = "L/D Aqua - Viability",
                              gateData = c(0, 0, 250000, 250000,
                                           0, 650, 650, 0)
                            )
                          )
      )

pipL <-
  addProcessingStep(
    pipL,
    whichQueue = "pre-processing",
    CytoProcessingStep(
      name = "perform_QC",
      FUN = "qualityControlPeacoQC",
      ARGS = list(
        preTransform = TRUE,
        min_cells = 150, # default
        max_bins = 500, # default
        step = 500, # default,
        MAD = 6, # default
        IT_limit = 0.55, # default
        force_IT = 150, # default
        peak_removal = 0.3333, # default
        min_nr_bins_peakdetection = 10 # default
      )
    )
  )

pipL <-
  addProcessingStep(pipL,
                    whichQueue = "pre-processing",
                    CytoProcessingStep(
                      name = "transform",
                      FUN = "applyScaleTransforms",
                      ARGS = list()
                    )
  )

# execute pipeline, remove cache if existing with the same experiment name
suppressWarnings(execute(pipL, rmCache = TRUE, path = outputDir))

# re-execute as is without removing cache => all results found in cache!
suppressWarnings(execute(pipL, rmCache = FALSE, path = outputDir))

### *** EXAMPLE 2: building CytoPipeline from JSON file *** ###

```

```

jsonDir <- system.file("extdata", package = "CytoPipeline")
jsonPath <- file.path(jsonDir, "pipelineParams.json")

pipl2 <- CytoPipeline(jsonPath,
                      experimentName = experimentName,
                      sampleFiles = sampleFiles)

# note we temporarily set working directory into package root directory
# needed as json path mentions "./" path for sample files
suppressWarnings(execute(pipl2, rmCache = TRUE, path = outputDir))

### *** EXAMPLE 3: building CytoPipeline from cache (previously run) *** ###

experimentName <- "OMIP021_PeacoQC"
pipl3 <- buildCytoPipelineFromCache(
  experimentName = experimentName,
  path = outputDir)

suppressWarnings(execute(pipl3,
  rmCache = FALSE,
  path = outputDir))

```

---

|                    |                                       |
|--------------------|---------------------------------------|
| exportCytoPipeline | <i>exporting CytoPipeline objects</i> |
|--------------------|---------------------------------------|

---

## Description

functions to export CytoPipeline objects in various formats

## Usage

```
export2JSONFile(x, path)
```

## Arguments

|      |                                                     |
|------|-----------------------------------------------------|
| x    | a CytoPipeline object                               |
| path | the full path to the name of the file to be created |

## Value

- for export2JSONFile: nothing

## Functions

- export2JSONFile(): exports a CytoPipeline object to a JSON file (writing the file = side effect)

## Examples

```
outputDir <- base::tempdir()

rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
experimentName <- "OMIP021_PeacoQC"
sampleFiles <- file.path(rawDataDir, list.files(rawDataDir,
                                                pattern = "Donor"))

# build CytoPipeline object using json input
jsonPath <- file.path(system.file("extdata", package = "CytoPipeline"),
                      "pipelineParams.json")

pipL <- CytoPipeline(jsonPath,
                    experimentName = experimentName,
                    sampleFiles = sampleFiles)

# remove the last pre-processing step
nPreProcessing <- getNbProcessingSteps(pipL, whichQueue = "pre-processing")
pipL <- removeProcessingStep(pipL, whichQueue = "pre-processing",
                           index = nPreProcessing)

# export back to json file
export2JSONFile(pipL, path = file.path(outputDir, "newFile.json"))
```

---

|                 |                                               |
|-----------------|-----------------------------------------------|
| findTimeChannel | <i>find time channel in flowSet/flowFrame</i> |
|-----------------|-----------------------------------------------|

---

## Description

tries to find a channel in a flowSet/flowFrame that could be the time channel. First tries to identify a channel name containing the 'time' string, then tries to identify a single monotonically increasing channel.

## Usage

```
findTimeChannel(obj, excludeChannels = c())
```

## Arguments

obj                    a flowCore::flowFrame or flowCore::flowSet  
 excludeChannels       vector of column names to exclude in the search

## Value

a character, name of the found channel that should be representing time. If not found, returns NULL.

**Examples**

```
data(OMIP021Samples)

ret <- findTimeChannel(OMIP021Samples[[1]])
ret # "Time"
```

---

```
getAcquiredCompensationMatrix
```

*extract compensation matrix from a flowCore::flowFrame*

---

**Description**

helper function retrieving the compensation matrix stored in fcs file (if any). It scans the following keywords: \$SPILL, \$spillover and \$SPILLOVER

**Usage**

```
getAcquiredCompensationMatrix(ff)
```

**Arguments**

ff                    a flowCore::flowFrame

**Value**

the found compensation matrix

**Examples**

```
rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
sampleFiles <-
  file.path(rawDataDir, list.files(rawDataDir, pattern = "Donor"))

truncateMaxRange <- FALSE
minLimit <- NULL

# create flowCore::flowSet with all samples of a dataset
fsRaw <- readSampleFiles(
  sampleFiles = sampleFiles,
  whichSamples = "all",
  truncate_max_range = truncateMaxRange,
  min.limit = minLimit)
compensationMatrix <- getAcquiredCompensationMatrix(fsRaw[[2]])
```

---

```
getChannelNamesFromMarkers
```

*get channel names from markers*

---

## Description

finds name of channels corresponding to user provided markers

## Usage

```
getChannelNamesFromMarkers(ff, markers)
```

## Arguments

|         |                                                                                                                                                                                                                                                                                      |
|---------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ff      | a flowCore::flowFrame                                                                                                                                                                                                                                                                |
| markers | a vector of markers, either provided as : <ul style="list-style-type: none"><li>• an array of booleans (referring to flowFrame columns)</li><li>• an array of integers (indices in flowFrame columns)</li><li>• an array of characters (exact markers or channel patterns)</li></ul> |

## Value

a character vector, containing the names of the corresponding channels

## Examples

```
data(OMIP021Samples)

# with existing markers
ret <- getChannelNamesFromMarkers(
  OMIP021Samples[[1]],
  c(
    "FSC-A",
    "L/D Aqua - Viability",
    "FITC - gdTCR",
    "PECy5 - CD28"
  )
)

ret # c("FSC-A", "525/50Violet-A", "530/30Blue-A", "670/30Yellow-A")

# with boolean vector
indices <- c(1, 6, 14, 18)
boolInput <- rep(FALSE, 21)
boolInput[indices] <- TRUE
ret2 <- getChannelNamesFromMarkers(
  OMIP021Samples[[1]],
  boolInput
)
```

```
ret2 # c("FSC-A", "525/50Violet-A", "530/30Blue-A", "670/30Yellow-A")

# with indices vector
ret3 <- getChannelNamesFromMarkers(
  OMIP021Samples[[1]],
  indices
)
ret3 # c("FSC-A", "525/50Violet-A", "530/30Blue-A", "670/30Yellow-A")
```

---

|                |                          |
|----------------|--------------------------|
| getFCSFileName | <i>get fcs file name</i> |
|----------------|--------------------------|

---

### Description

get basename of \$FILENAME keyword if exists

### Usage

```
getFCSFileName(ff)
```

### Arguments

ff                      a flowCore::flowFrame

### Value

the basename of \$FILENAME keyword

### Examples

```
data(OMIP021Samples)

fName <- getFCSFileName(OMIP021Samples[[1]])
```

---

|                  |                                                             |
|------------------|-------------------------------------------------------------|
| getTransfoParams | <i>get transformation parameters for a specific channel</i> |
|------------------|-------------------------------------------------------------|

---

### Description

investigates a flowCore::transformList object to get the type and parameters of the transformation applying to a specific channel

### Usage

```
getTransfoParams(transList, channel)
```

**Arguments**

|           |                           |
|-----------|---------------------------|
| translist | a flowCore::transformList |
| channel   | channel name              |

**Value**

If the transformation exists for the specified channel, and is either recognized as a logicle transfo or a linear transfo, a list with two slots:

- \$type a character containing the transfo type ('logicle' or 'linear')
- \$params\_list a list of named numeric, according to transfo type

Otherwise, NULL is returned.

**Examples**

```
data(OMIP021Samples)

# set-up a hybrid transformation list :
# - two channels are logicle-ly transformed with automatic param estimates
# - one channel has explicit logicle transfo with default parameters
# - one channel has linear transformation
# - other channels have no transformation
translist <- flowCore::estimateLogicle(
  OMIP021Samples[[1]],
  c("450/50Violet-A", "525/50Violet-A")
)
translist <- c(
  translist,
  flowCore::transformList(
    "FSC-A",
    flowCore::linearTransform(
      a = 0.1,
      b = 0
    )
  ),
  flowCore::transformList(
    "540/30Violet-A",
    flowCore::logicleTransform()
  )
)

ret1 <- getTransfoParams(translist, channel = "FSC-A")
ret1$type # "linear"
ret1$paramsList # a = 0.1, b = 0.

ret2 <- getTransfoParams(translist, channel = "525/50Violet-A")
ret2$type # "logicle"
ret2$paramsList # a = 0., w = 0.2834, m = 4.5, t = 262143

ret3 <- getTransfoParams(translist, channel = "540/30Violet-A")
```

```
ret3$type # "logicle"
ret3$paramsList # a = 0., w = 0.5, m = 4.5, t = 262144
```

---

|              |                                               |
|--------------|-----------------------------------------------|
| ggplotEvents | <i>plot events in 1D or 2D, using ggplot2</i> |
|--------------|-----------------------------------------------|

---

## Description

plot events of specific channels of either : flowCore::flowFrame, or flowCore::flowSet in 2D or 1D, mimicking FlowJo type of graph.  
 if 1D : geom\_density will be used  
 if 2D : geom\_hex will be used

## Usage

```
ggplotEvents(
  obj,
  xChannel,
  yChannel = NULL,
  nDisplayCells = Inf,
  seed = NULL,
  bins = 216,
  fill = "lightblue",
  alpha = 0.2,
  xScale = c("linear", "logicle"),
  yScale = c("linear", "logicle"),
  xLogicleParams = NULL,
  yLogicleParams = NULL,
  xLinearRange = NULL,
  yLinearRange = NULL,
  transList = NULL,
  runTransforms = FALSE
)
```

## Arguments

|               |                                                                                                                             |
|---------------|-----------------------------------------------------------------------------------------------------------------------------|
| obj           | a flowCore::flowFrame or flowCore::flowSet                                                                                  |
| xChannel      | channel (name or index) or marker name to be displayed on x axis                                                            |
| yChannel      | channel (name or index) or marker name to be displayed on y axis                                                            |
| nDisplayCells | maximum number of events that will be plotted. If the number of events exceed this number, a sub-sampling will be performed |
| seed          | seed used for sub-sampling (if any)                                                                                         |
| bins          | used in geom_hex                                                                                                            |
| fill          | used in geom_density                                                                                                        |

|                |                                                                                                                                                                                                                                                                                                                                    |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| alpha          | used in geom_density                                                                                                                                                                                                                                                                                                               |
| xScale         | scale to be used for the x axis (note "linear" corresponds to no transformation)                                                                                                                                                                                                                                                   |
| yScale         | scale to be used for the y axis (note "linear" corresponds to no transformation)                                                                                                                                                                                                                                                   |
| xLogicleParams | if (xScale == "logicle"), the parameters of the logicle transformation to be used, as a list(w = ..., m = ..., a = ..., t = ...). If NULL, these parameters will be estimated by flowCore::estimateLogicle()                                                                                                                       |
| yLogicleParams | if (yScale == "logicle"), the parameters of the logicle transformation to be used, as a list(w = ..., m = ..., a = ..., t = ...). If NULL, these parameters will be estimated by flowCore::estimateLogicle()                                                                                                                       |
| xLinearRange   | if (xScale == "linear"), the x axis range to be used                                                                                                                                                                                                                                                                               |
| yLinearRange   | if (yScale == "linear"), the y axis range to be used                                                                                                                                                                                                                                                                               |
| transList      | optional list of scale transformations to be applied to each channel. If it is non null, 'x/yScale', 'x/yLogicleParams' and 'x/yLinear_range' will be discarded.                                                                                                                                                                   |
| runTransforms  | (TRUE/FALSE) Will the application of non linear scale result in data being effectively transformed ? <ul style="list-style-type: none"> <li>• If TRUE, than the data will undergo transformations prior to visualization.</li> <li>• If FALSE, the axis will be scaled but the data themselves will not be transformed.</li> </ul> |

## Value

a list of ggplot objects

## Examples

```
data(OMIP021Samples)

### 1D Examples ###

# simple linear scale example
ggplotEvents(OMIP021Samples[[1]],
             xChannel = "FSC-A",
             xScale = "linear")

# with explicit linear range
ggplotEvents(OMIP021Samples[[1]],
             xChannel = "FSC-A",
             xScale = "linear",
             xLinearRange = c(0, 250000))

# with linear scale, several flow frames
ggplotEvents(OMIP021Samples, xChannel = "FSC-A", xScale = "linear")

# simple logicle scale example
ggplotEvents(OMIP021Samples[[1]],
             xChannel = "450/50Violet-A",
             xScale = "logicle")
```

```

# logicle scale, explicit parameters
ggplotEvents(OMIP021Samples[[1]],
              xChannel = "450/50Violet-A",
              xScale = "logicle", xLogicleParams = list(
                a = 1,
                w = 2,
                m = 7,
                t = 270000))

# with sub-sampling
ggplotEvents(OMIP021Samples[[2]],
              xChannel = "450/50Violet-A",
              xScale = "logicle", nDisplayCells = 5000)

# tuning some plot parameters
ggplotEvents(OMIP021Samples[[2]],
              xChannel = "450/50Violet-A",
              xScale = "logicle", alpha = 0.5, fill = "red")

# examples that use a transformation list, estimated after compensation
compensationMatrix <- flowCore::spillover(OMIP021Samples[[1]])$SPILL

ffc <- runCompensation(OMIP021Samples[[1]],
                      spillover = compensationMatrix,
                      updateChannelNames = FALSE)

transList <- flowCore::estimateLogicle(
  ffc,
  colnames(compensationMatrix))

transList <-
  c(transList,
    flowCore::transformList(
      "FSC-A",
      flowCore::linearTransform(a = 0.00001)))

# linear example, without running the transformations on data
ggplotEvents(OMIP021Samples[[1]],
              xChannel = "450/50Violet-A",
              xScale = "linear",
              transList = transList,
              runTransforms = FALSE)

# linear example, now running the transformations on data
ggplotEvents(OMIP021Samples[[1]],
              xChannel = "450/50Violet-A",
              xScale = "linear",
              transList = transList,
              runTransforms = TRUE)

# logicle example, without running the transformations on data
ggplotEvents(OMIP021Samples[[1]],
              xChannel = "FSC-A",

```

```
xScale = "logicle",
transList = transList,
runTransforms = FALSE)

# logicle example, now running the transformations on data
ggplotEvents(OMIP021Samples[[1]],
  xChannel = "FSC-A",
  xScale = "logicle",
  transList = transList,
  runTransforms = TRUE)

### 2D examples ###

# simple linear example
ggplotEvents(OMIP021Samples[[1]],
  xChannel = "FSC-A",
  xScale = "linear",
  yChannel = "610/20Violet-A",
  yScale = "logicle")

# simple linear example, 2 flow frames
ggplotEvents(OMIP021Samples,
  xChannel = "FSC-A",
  xScale = "linear",
  yChannel = "SSC-A",
  yScale = "linear")

# logicle vs linear example
ggplotEvents(OMIP021Samples[[1]],
  xChannel = "450/50Violet-A",
  xScale = "logicle",
  yChannel = "SSC-A",
  yScale = "linear")

# 2X logicle example
ggplotEvents(OMIP021Samples[[1]],
  xChannel = "TETaGC",
  xScale = "logicle",
  yChannel = "CD27",
  yScale = "logicle")

# tuning nb of bins
ggplotEvents(OMIP021Samples[[1]],
  xChannel = "TETaGC",
  xScale = "logicle",
  yChannel = "CD27",
  yScale = "logicle",
  bins = 128)

# using transformation list, not run on data
ggplotEvents(OMIP021Samples[[1]],
  xChannel = "TETaGC",
```

```

        xScale = "logicle",
        yChannel = "CD27",
        yScale = "logicle",
        transList = transList,
        runTransforms = FALSE)

# using transformation list, run on data
ggplotEvents(OMIP021Samples[[1]],
             xChannel = "TETaGC",
             xScale = "logicle",
             yChannel = "CD27",
             yScale = "logicle",
             transList = transList,
             runTransforms = TRUE)

```

---

ggplotFilterEvents      *plot filtered events in 2D, using ggplot*

---

## Description

plot events of specific channels of either : flowCore::flowFrame, or flowCore::flowSet in 2D, showing the impact of applying a filter between :

- a 'pre' flowframe

## Usage

```

ggplotFilterEvents(
  ffPre,
  ffPost,
  xChannel,
  yChannel,
  nDisplayCells = 10000,
  seed = NULL,
  size = 0.5,
  xScale = c("linear", "logicle"),
  yScale = c("linear", "logicle"),
  xLogicleParams = NULL,
  yLogicleParams = NULL,
  xLinearRange = NULL,
  yLinearRange = NULL,
  transList = NULL,
  runTransforms = FALSE,
  interactive = FALSE
)

```

**Arguments**

|                             |                                                                                                                                                                                                                                                                                                                                |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <code>ffPre</code>          | a <code>flowCore::flowFrame</code> , before applying filter                                                                                                                                                                                                                                                                    |
| <code>ffPost</code>         | a <code>flowCore::flowFrame</code> , after applying filter                                                                                                                                                                                                                                                                     |
| <code>xChannel</code>       | channel (name or index) or marker name to be displayed on x axis                                                                                                                                                                                                                                                               |
| <code>yChannel</code>       | channel (name or index) or marker name to be displayed on y axis                                                                                                                                                                                                                                                               |
| <code>nDisplayCells</code>  | maximum number of events that will be plotted. If the number of events exceed this number, a subsampling will be performed                                                                                                                                                                                                     |
| <code>seed</code>           | seed used for sub-sampling (if any)                                                                                                                                                                                                                                                                                            |
| <code>size</code>           | used by <code>geom_point()</code>                                                                                                                                                                                                                                                                                              |
| <code>xScale</code>         | scale to be used for the x axis (note "linear" corresponds to no transformation)                                                                                                                                                                                                                                               |
| <code>yScale</code>         | scale to be used for the y axis (note "linear" corresponds to no transformation)                                                                                                                                                                                                                                               |
| <code>xLogicleParams</code> | if ( <code>xScale == "logicle"</code> ), the parameters of the logicle transformation to be used, as a list( <code>w = ..., m = ..., a = ..., t = ...</code> ) If <code>NULL</code> , these parameters will be estimated by <code>flowCore::estimateLogicle()</code>                                                           |
| <code>yLogicleParams</code> | if ( <code>yScale == "logicle"</code> ), the parameters of the logicle transformation to be used, as a list( <code>w = ..., m = ..., a = ..., t = ...</code> ) If <code>NULL</code> , these parameters will be estimated by <code>flowCore::estimateLogicle()</code>                                                           |
| <code>xLinearRange</code>   | if ( <code>xScale == "linear"</code> ), linear range to be used                                                                                                                                                                                                                                                                |
| <code>yLinearRange</code>   | if ( <code>yScale == "linear"</code> ), linear range to be used                                                                                                                                                                                                                                                                |
| <code>transList</code>      | optional list of scale transformations to be applied to each channel. If it is non null, <code>'x/yScale'</code> , <code>'x/yLogicleParams'</code> and <code>'x/yLinear_range'</code> will be discarded.                                                                                                                       |
| <code>runTransforms</code>  | (TRUE/FALSE) Will the application of non linear scale result in data being effectively transformed ? <ul style="list-style-type: none"> <li>• If TRUE, than the data will undergo transformations prior to visualization.</li> <li>• If FALSE, the axis will be scaled but the data themselves are not transformed.</li> </ul> |
| <code>interactive</code>    | if TRUE, transform the scaling formats such that the <code>ggcyto::x_scale_logicle()</code> and <code>ggcyto::y_scale_logicle()</code> do work with <code>plotly::ggplotly()</code>                                                                                                                                            |

**Value**

a ggplot object

**Examples**

```
data(OMIP021Samples)

ffPre <- OMIP021Samples[[1]]

# creating a manual polygon gate filter based on channels L/D and FSC-A

LDMarker <- "L/D Aqua - Viability"

LDChannel <- getChannelNamesFromMarkers(ffPre, markers = LDMarker)
```

```

liveGateMatrix <- matrix(
  data = c(
    50000, 50000, 100000, 200000, 200000,
    100, 1000, 2000, 2000, 1
  ),
  ncol = 2,
  dimnames = list(
    c(),
    c("FSC-A", LDChannel)
  )
)

liveGate <- flowCore::polygonGate(
  filterId = "Live",
  .gate = liveGateMatrix
)

selectedLive <- flowCore::filter(ffPre, liveGate)
ffL <- flowCore::Subset(ffPre, selectedLive)

# show the results

# subsample 5000 points
ggplotFilterEvents(
  ffPre = ffPre,
  ffPost = ffL,
  nDisplayCells = 5000,
  xChannel = "FSC-A", xScale = "linear",
  yChannel = LDMarker, yScale = "logicle") +
  ggplot2::ggtitle("Live gate filter - 5000 points")

# with all points
ggplotFilterEvents(
  ffPre = ffPre,
  ffPost = ffL,
  nDisplayCells = Inf,
  xChannel = "FSC-A", xScale = "linear",
  yChannel = LDMarker, yScale = "logicle") +
  ggplot2::ggtitle("Live gate filter - all points")

```

---

ggplotFlowRate

*plot flow rate as a function of time, using ggplot2*


---

## Description

plot flow rate as a function of time, using ggplot2

**Usage**

```
ggplotFlowRate(obj, title = "Flow Rate", timeUnit = 100)
```

**Arguments**

|          |                                                                                 |
|----------|---------------------------------------------------------------------------------|
| obj      | a flowCore::flowFrame or flowCore::flowSet                                      |
| title    | a title for the graph                                                           |
| timeUnit | which time interval is used to calculate "instant" flow rate (default = 100 ms) |

**Value**

a ggplot graph

**Examples**

```
data(OMIP021Samples)

# single flowFrame plot
ggplotFlowRate(OMIP021Samples[[1]])

# two flowFrames plot
ggplotFlowRate(OMIP021Samples)

# single plot with title
ggplotFlowRate(OMIP021Samples[[1]], title = "Test Flow Rate plot")

# explicit time unit
ggplotFlowRate(OMIP021Samples[[1]], timeUnit = 50)
```

---

handlingProcessingSteps

*handling processing steps in CytoPipeline objects*

---

**Description**

functions to manipulate processing steps in processing queues of CytoPipeline objects

**Usage**

```
addProcessingStep(
  x,
  whichQueue = c("scale transform", "pre-processing"),
  newPS
)

removeProcessingStep(
```

```

    x,
    whichQueue = c("scale transform", "pre-processing"),
    index
)

getNbProcessingSteps(x, whichQueue = c("scale transform", "pre-processing"))

getProcessingStep(
  x,
  whichQueue = c("scale transform", "pre-processing"),
  index
)

getProcessingStepNames(x, whichQueue = c("scale transform", "pre-processing"))

cleanProcessingSteps(
  x,
  whichQueue = c("both", "scale transform", "pre-processing")
)

showProcessingSteps(x, whichQueue = c("scale transform", "pre-processing"))

```

### Arguments

|                         |                                                                       |
|-------------------------|-----------------------------------------------------------------------|
| <code>x</code>          | a CytoPipeline object                                                 |
| <code>whichQueue</code> | selects the processing queue for which we manage the processing steps |
| <code>newPS</code>      | the new processing step to be added (CytoProcessingStep object)       |
| <code>index</code>      | index of the processing step to remove                                |

### Value

- for `addProcessingStep`: the updated CytoPipeline object
- for `removeProcessingStep`: the updated CytoPipeline object
- for `getNbProcessingSteps`: the number of processing steps present in the target queue
- for `getProcessingStep`: the obtained CytoProcessingStep object
- for `getProcessingStepNames`: the vector of step names
- for `cleanProcessingSteps`: the updated CytoPipeline object
- for `showProcessingSteps`: nothing (only console display side effect is required)

## Functions

- `addProcessingStep()`: adds a processing step in one of the processing queues (at the end), returns the modified `CytoPipeline` object
- `removeProcessingStep()`: removes a processing step from one of the processing queues, returns the modified `CytoPipeline` object
- `getNbProcessingSteps()`: gets the number of processing steps in a processing queue
- `getProcessingStep()`: gets a processing step at a specific index of a processing queue
- `getProcessingStepNames()`: gets a character vector of all processing step names of a specific processing queue
- `cleanProcessingSteps()`: deletes all processing steps in one or both processing queues, returns the modified `CytoPipeline` object
- `showProcessingSteps()`: shows all processing steps in a processing queue

## Examples

```
rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
experimentName <- "OMIP021_PeacoQC"
sampleFiles <- file.path(rawDataDir, list.files(rawDataDir,
                                                pattern = "Donor"))

transListPath <-
  file.path(system.file("extdata", package = "CytoPipeline"),
            "OMIP021_TransList.rds")

# main parameters : sample files and experiment name
pipelineParams <- list()
pipelineParams$experimentName <- experimentName
pipelineParams$sampleFiles <- sampleFiles

# create CytoPipeline object (no step defined yet)
pipL <- CytoPipeline(pipelineParams)

# add a processing step in scale transformation queue
pipL <- addProcessingStep(pipL,
  whichQueue = "scale transform",
  CytoProcessingStep(
    name = "scale_transform_read",
    FUN = "readRDS",
    ARGS = list(file = transListPath)
  ))

getNbProcessingSteps(pipL, "scale transform") # returns 1

# add another processing step in scale transformation queue
pipL <- addProcessingStep(pipL,
  whichQueue = "scale transform",
  CytoProcessingStep(
    name = "scale_transform_sum",
    FUN = "sum",
```

```

        ARGs = list()
    )

    getNbProcessingSteps(pipL, "scale transform") # returns 2

    getProcessingStepNames(pipL, whichQueue = "scale transform")

    # removes second processing step in scale transformation queue
    pipL <- removeProcessingStep(pipL,
                                whichQueue = "scale transform",
                                index = 2)

    # get processing step object
    pS <- getProcessingStep(pipL, whichQueue = "scale transform", index = 1)
    getCPSName(pS) #"scale_transform_read"

    # add a processing step in pre-processing queue
    pipL <- addProcessingStep(pipL,
                              whichQueue = "pre-processing",
                              CytoProcessingStep(
                                  name = "pre-processing_sum",
                                  FUN = "sum",
                                  ARGs = list()
                              ))

    getNbProcessingSteps(pipL, "scale transform") # returns 1
    getNbProcessingSteps(pipL, "pre-processing") # returns also 1

    showProcessingSteps(pipL, whichQueue = "scale transform")
    showProcessingSteps(pipL, whichQueue = "pre-processing")

    # cleans both processing queues
    pipL <- cleanProcessingSteps(pipL)
    pipL

```

---

inspectCytoPipelineObjects

*inspect CytoPipeline results objects*

---

## Description

functions to obtain results objects formats

## Usage

```

getCytoPipelineExperimentNames(
  path = ".",
  pattern = NULL,
  ignore.case = FALSE,
  fixed = FALSE

```

```
)

getCytoPipelineObjectFromCache(
  x,
  path = ".",
  whichQueue = c("scale transform", "pre-processing"),
  sampleFile = NULL,
  objectName
)

getCytoPipelineObjectInfos(
  x,
  path = ".",
  whichQueue = c("scale transform", "pre-processing"),
  sampleFile = NULL
)

getCytoPipelineFlowFrame(
  x,
  path = ".",
  whichQueue = c("scale transform", "pre-processing"),
  sampleFile,
  objectName
)

getCytoPipelineScaleTransform(
  x,
  path = ".",
  whichQueue = c("scale transform", "pre-processing"),
  sampleFile = NULL,
  objectName
)

plotCytoPipelineProcessingQueue(
  x,
  whichQueue = c("pre-processing", "scale transform"),
  purpose = c("run status", "description"),
  sampleFile = NULL,
  path = ".",
  title = TRUE,
  box.type = "ellipse",
  lwd = 1,
  box.prop = 0.5,
  box.cex = 0.7,
  cex.txt = 0.7,
  box.size = 0.1,
  dtext = 0.15,
  ...
)
```

)

```
collectNbOfRetainedEvents(experimentName, path = ".", whichSampleFiles)
```

### Arguments

|                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| path             | root path to locate the search for file caches                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| pattern          | optional pattern limiting the search for experiment names                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| ignore.case      | (TRUE/FALSE) used in pattern matching (grepl)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| fixed            | (TRUE/FALSE) used in pattern matching (grepl)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| x                | a CytoPipeline object                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| whichQueue       | which queue to look into                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| sampleFile       | which sampleFile is looked for: <ul style="list-style-type: none"> <li>• if whichQueue == "scale transform", the sampleFile is ignored</li> <li>• if NULL and whichQueue == "pre-processing", the sampleFile is defaulted to the first one belonging to the experiment</li> </ul>                                                                                                                                                                                                                                                                                                                |
| objectName       | (character) which object name to look for                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| purpose          | purpose of the workflow plot <ul style="list-style-type: none"> <li>• if "run status" (default), the disk cache will be inspected and the box colours will be set according to run status (green = run, orange = not run, red = definition not consistent with cache). Moreover, the object classes and names will be filled in if found in the cache.</li> <li>• if "description", the workflow will be obtained from the step definition in the x object, not from the disk cache. As a result, all boxes will be coloured in black, and no object class and name will be provided.</li> </ul> |
| title            | if TRUE, adds a title to the plot                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| box.type         | shape of label box (rect, ellipse, diamond, round, hexa, multi)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| lwd              | default line width of arrow and box (one numeric value)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| box.prop         | length/width ratio of label box (one numeric value)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| box.cex          | relative size of text in boxes (one numeric value)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| cex.txt          | relative size of arrow text (one numeric value)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| box.size         | size of label box (one numeric value)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| dtext            | controls the position of arrow text relative to arrowhead (one numeric value)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| ...              | other arguments passed to diagram::plotmat()                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| experimentName   | the experimentName used to select the file cache on disk                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| whichSampleFiles | indicates for which sample files the number of retained events are to be collected. If missing, all sample files will be used.                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |

**Value**

- for `getCytoPipelineExperimentNames`: a vector of character containing found experiment names
- for `getCytoPipelineObjectFromCache`: the found object (or stops with an error message if the target object is not found)
- for `getCytoPipelineObjectInfos`: a dataframe with the collected information about the found objects (or stops with an error message if no target object was found)
- for `getCytoPipelineFlowFrame`: the found flowFrame (or stops with an error message if the target object is not found, or if the object is no flowFrame)
- for `getCytoPipelineScaleTransform`: the found flowFrame (or stops with an error message if the target object is not found, or if the object is no transformList)
- for `plotCytoPipelineProcessingQueue`: nothing
- for `collectNbOfRetainedEvents`: a dataframe with the collected number of events columns refer to pre-processing steps rows refer to samples

**Functions**

- `getCytoPipelineExperimentNames()`: This function looks into a path for stored file caches and gets the corresponding experiment names
- `getCytoPipelineObjectFromCache()`: Given a CytoPipeline object, this function retrieves a specific object in the corresponding file cache
- `getCytoPipelineObjectInfos()`: Given a CytoPipeline object, this function retrieves the information related to a specific object name, i.e. object name and object class
- `getCytoPipelineFlowFrame()`: Given a CytoPipeline object, this function retrieves a specific flowCore::flowFrame object in the corresponding file cache object name and object class
- `getCytoPipelineScaleTransform()`: Given a CytoPipeline object, this function retrieves a specific flowCore::transformList object in the corresponding file cache
- `plotCytoPipelineProcessingQueue()`: This functions displays a plot of a processing queue of a CytoPipeline object, using `diagram::plotmat()`.
  - If a step is in run state for all sample files, the corresponding box appears in green
  - If a step is in non run state for at least one sample file, the corresponding box appears in orange
  - If at least one step is not consistent with cache, the whole set of boxes appears in red
- `collectNbOfRetainedEvents()`: Given a CytoPipeline object, this function retrieves, for all pre-processing steps, given the output is a flowFrame, the number of retained event.

[illegible]

```

        path = outputDir,
        sampleFile = sampleFiles(pipL2)[1])

# get any object (output of one step)
obj <-
  getCytoPipelineObjectFromCache(pipL2, whichQueue = "scale transform",
                                objectName = "compensate_obj",
                                path = outputDir)

class(obj) # flowCore::flowSet

# collect number of retained events at each step
nbEventsDF <- collectNbOfRetainedEvents(
  experimentName = experimentNames[1],
  path = outputDir)

```

---

interactingWithCytoPipelineCache

*interaction between CytoPipeline object and disk cache*


---

## Description

functions supporting the interaction between a CytoPipeline object and the file cache on disk

## Usage

```

deleteCytoPipelineCache(x, path = ".")

buildCytoPipelineFromCache(experimentName, path = ".")

checkCytoPipelineConsistencyWithCache(
  x,
  path = ".",
  whichQueue = c("both", "scale transform", "pre-processing"),
  sampleFile = NULL
)

```

## Arguments

|                |                                                                                              |
|----------------|----------------------------------------------------------------------------------------------|
| x              | a CytoPipeline object                                                                        |
| path           | the full path to the experiment storage on disk (without the /.cache)                        |
| experimentName | the experimentName used to select the file cache on disk                                     |
| whichQueue     | which processing queue to check the consistency of                                           |
| sampleFile     | if whichQueue == "pre-processing" or "both": which sample file(s) to check on the disk cache |

**Value**

for deleteCytoPipelineCache: TRUE if successfully removed  
 for buildCytoPipelineFromCache: the built CytoPipeline object  
 for checkCytoPipelineConsistencyWithCache: a list with the following values:

- isConsistent (TRUE/FALSE)
- inconsistencyMsg: character filled in by an inconsistency message in case the cache and CytoPipeline object are not consistent with each other
- scaleTransformStepStatus: a character vector, containing, for each scale transform step, a status from c("run", "not run", "inconsistent")
- preProcessingStepStatus: a character matrix, containing, for each pre-processing step (rows), for each sample file (columns), a status from c("run", "not run", "inconsistent")

**Functions**

- deleteCytoPipelineCache(): delete the whole disk cache corresponding to the experiment of a CytoPipeline object
- buildCytoPipelineFromCache(): builds a new CytoPipeline object, based on the information stored in the file cache
- checkCytoPipelineConsistencyWithCache(): check the consistency between the processing steps described in a CytoPipeline object, and what is stored in the file cache

**Examples**

```
# preliminary run:
# build CytoPipeline object using json input, run and store results in cache
rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
experimentName <- "OMIP021_PeacoQC"
sampleFiles <- file.path(rawDataDir, list.files(rawDataDir,
  pattern = "Donor"))

jsonDir <- system.file("extdata", package = "CytoPipeline")
jsonPath <- file.path(jsonDir, "pipelineParams.json")
outputDir <- base::tempdir()
pipL <- CytoPipeline(jsonPath,
  experimentName = experimentName,
  sampleFiles = sampleFiles)

# note we temporarily set working directory into package root directory
# needed as json path mentions "../" path for sample files
suppressWarnings(execute(pipL, rmCache = TRUE, path = outputDir))

# rebuild CytoPipeline from stored results in cache, for a specific
# experiment

experimentName <- "OMIP021_PeacoQC"
pipL2 <- buildCytoPipelineFromCache(
```

```

    experimentName = experimentName,
    path = outputDir)

# checking consistency between CytoPipeline object and cache
res <- checkCytoPipelineConsistencyWithCache(pipL2)
#res

suppressWarnings(execute(pipL2, rmCache = FALSE, path = outputDir))
# (everything is already stored in cache)

# deleting cache related to a specific experiment
pipL3 <- CytoPipeline(experimentName = experimentName)
deleteCytoPipelineCache(pipL3, path = outputDir)

```

---

OMIP021Samples

*OMIP021Samples dataset*


---

### Description

OMIP021Samples dataset

### Format

a flowCore::flowSet with two different flowFrames each one contains one flow cytometry sample corresponding to Donor1.fcs and Donor2.fcs in following source. A subsampling of 5,000 events has been performed on each file.

### Value

nothing

### Source

<https://flowrepository.org/experiments/305>

---

qualityControlFlowAI

*perform QC with flowAI*


---

### Description

this function is a wrapper around flowAI::flow\_auto\_qc() function. It also pre-selects the channels to be handled (=> all signal channels)

**Usage**

```
qualityControlFlowAI(
  ff,
  preTransform = FALSE,
  transList = NULL,
  outputDiagnostic = FALSE,
  outputDir = NULL,
  ...
)
```

**Arguments**

|                               |                                                                                                       |
|-------------------------------|-------------------------------------------------------------------------------------------------------|
| <code>ff</code>               | a <code>flowCore::flowFrame</code>                                                                    |
| <code>preTransform</code>     | if TRUE, apply the <code>transList</code> scale transform prior to running the gating algorithm       |
| <code>transList</code>        | applied in conjunction with <code>preTransform</code>                                                 |
| <code>outputDiagnostic</code> | if TRUE, stores diagnostic files generated by <code>flowAI</code> in <code>outputDir</code> directory |
| <code>outputDir</code>        | used in conjunction with <code>outputDiagnostic</code>                                                |
| <code>...</code>              | additional parameters passed to <code>flowAI::flow_auto_qc()</code>                                   |

**Value**

a `flowCore::flowFrame` with removed low quality events from the input

**Examples**

```
rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
sampleFiles <-
  file.path(rawDataDir, list.files(rawDataDir, pattern = "Donor"))

truncateMaxRange <- FALSE
minLimit <- NULL

# create flowCore::flowSet with all samples of a dataset
fsRaw <- readSampleFiles(
  sampleFiles = sampleFiles,
  whichSamples = "all",
  truncate_max_range = truncateMaxRange,
  min.limit = minLimit)

suppressWarnings(ff_QualityControl <-
  qualityControlFlowAI(fsRaw[[2]],
    remove_from = "all", # all default
    second_fractionFR = 0.1,
    deviationFR = "MAD",
    alphaFR = 0.01,
    decompFR = TRUE,
```

```

outlier_binsFS = FALSE,
pen_valueFS = 500,
max_cptFS = 3,
sideFM = "both",
neg_valuesFM = 1))

```

---

qualityControlPeacoQC *perform QC with PeacoQC*

---

### Description

this function is a wrapper around `PeacoQC::PeacoQC()` function. It also pre-selects the channels to be handled (=> all signal channels)

### Usage

```

qualityControlPeacoQC(
  ff,
  preTransform = FALSE,
  transList = NULL,
  outputDiagnostic = FALSE,
  outputDir = NULL,
  ...
)

```

### Arguments

|                               |                                                                                                        |
|-------------------------------|--------------------------------------------------------------------------------------------------------|
| <code>ff</code>               | a <code>flowCore::flowFrame</code>                                                                     |
| <code>preTransform</code>     | if TRUE, apply the <code>transList</code> scale transform prior to running the gating algorithm        |
| <code>transList</code>        | applied in conjunction with <code>preTransform</code>                                                  |
| <code>outputDiagnostic</code> | if TRUE, stores diagnostic files generated by <code>PeacoQC</code> in <code>outputDir</code> directory |
| <code>outputDir</code>        | used in conjunction with <code>outputDiagnostic</code>                                                 |
| <code>...</code>              | additional parameters passed to <code>PeacoQC::PeacoQC()</code>                                        |

### Value

a `flowCore::flowFrame` with removed low quality events from the input

**Examples**

```

rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
sampleFiles <-
  file.path(rawDataDir, list.files(rawDataDir, pattern = "Donor"))

truncateMaxRange <- FALSE
minLimit <- NULL

# create flowCore::flowSet with all samples of a dataset
fsRaw <- readSampleFiles(
  sampleFiles = sampleFiles,
  whichSamples = "all",
  truncate_max_range = truncateMaxRange,
  min.limit = minLimit)

suppressWarnings(ff_m <- removeMarginsPeacoQC(x = fsRaw[[2]]))

ff_c <-
  compensateFromMatrix(ff_m,
    matrixSource = "fcs")

transList <-
  estimateScaleTransforms(
    ff = ff_c,
    fluoMethod = "estimateLogicle",
    scatterMethod = "linear",
    scatterRefMarker = "BV785 - CD3")

ff_QualityControl <- suppressWarnings(
  qualityControlPeacoQC(
    ff_c,
    preTransform = TRUE,
    transList = transList,
    min_cells = 150,
    max_bins = 500,
    MAD = 6,
    IT_limit = 0.55,
    force_IT = 150,
    peak_removal = (1/3),
    min_nr_bins_peakdetection = 10))

```

---

readRDSObject

*read RDS object*


---

**Description**

wrapper around readRDS, which discards any additional parameters passed in (...)

**Usage**

```
readRDSObject(RDSFile, ...)
```

**Arguments**

```
RDSFile      a RDS file containing a R object object
...          other arguments (not used)
```

**Value**

the read R object

**Examples**

```
data(OMIP021Samples)

transListPath <- file.path(system.file("extdata",
                                         package = "CytoPipeline"),
                           "OMIP021_TransList.rds")

transList <- readRDSObject(transListPath)

ff_c <- compensateFromMatrix(OMIP021Samples[[1]],
                             matrixSource = "fcs")

ff_t <- applyScaleTransforms(ff_c, transList = transList)
```

---

|                 |                              |
|-----------------|------------------------------|
| readSampleFiles | <i>Read fcs sample files</i> |
|-----------------|------------------------------|

---

**Description**

Wrapper around `flowCore::read.fcs()` or `flowCore::read.flowSet()`. Also adds a "Cell\_ID" additional column, used in `flowFrames` comparison

**Usage**

```
readSampleFiles(
  sampleFiles,
  whichSamples = "all",
  nSamples = NULL,
  seed = NULL,
  channelMarkerFile = NULL,
  ...
)
```

**Arguments**

|                                |                                                                                                                                                                                                                                                                                                                                                                                           |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <code>sampleFiles</code>       | a vector of character path to sample files                                                                                                                                                                                                                                                                                                                                                |
| <code>whichSamples</code>      | one of: <ul style="list-style-type: none"> <li>• 'all' if all sample files need to be read</li> <li>• 'random' if some samples need to be chosen randomly (in that case, using <code>nSamples</code> and <code>seed</code>)</li> <li>• a vector of indexes pointing to the <code>sampleFiles</code> vector</li> </ul>                                                                     |
| <code>nSamples</code>          | number of samples to randomly select (if <code>whichSamples == "random"</code> ). If <code>nSamples</code> is higher than nb of available samples, the output will be all samples                                                                                                                                                                                                         |
| <code>seed</code>              | an optional seed parameters (provided to ease reproducibility).                                                                                                                                                                                                                                                                                                                           |
| <code>channelMarkerFile</code> | an optional path to a csv file which provides the mapping between channels and markers. If provided, this csv file should contain a <code>Channel</code> column, and a <code>Marker</code> column. Optionally a 'Used' column can be provided as well (TRUE/FALSE). Channels for which the 'Used' column is set to FALSE will not be incorporated in the created <code>flowFrame</code> . |
| <code>...</code>               | additional parameters passed to <code>flowCore</code> file reading functions.                                                                                                                                                                                                                                                                                                             |

**Value**

either a `flowCore::flowSet` or a `flowCore::flowFrame` if `length(sampleFiles) == 1`

**Examples**

```
rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
sampleFiles <-
  file.path(rawDataDir, list.files(rawDataDir, pattern = "Donor"))

truncateMaxRange <- FALSE
minLimit <- NULL

# create flowCore::flowSet with all samples of a dataset
res <- readSampleFiles(
  sampleFiles = sampleFiles,
  whichSamples = "all",
  truncate_max_range = truncateMaxRange,
  min.limit = minLimit)

#res

# create a flowCore::flowFrame with one single sample
res2 <- readSampleFiles(
  sampleFiles = sampleFiles,
  whichSamples = 2,
  truncate_max_range = truncateMaxRange,
  min.limit = minLimit)
```

```
#res2
```

---

|                |                                         |
|----------------|-----------------------------------------|
| removeChannels | <i>remove channels from a flowFrame</i> |
|----------------|-----------------------------------------|

---

### Description

: in a flowCore::flowFrame, remove the channels of the given names.

### Usage

```
removeChannels(ff, channels)
```

### Arguments

|          |                                 |
|----------|---------------------------------|
| ff       | a flowCore::flowFrame           |
| channels | the channel names to be removed |

### Value

a new flowCore::flowFrame with the removed channels

### Examples

```
data(OMIP021Samples)

retFF <- removeChannels(OMIP021Samples[[1]],
  channel = "FSC-A")
```

---

|                           |                                                               |
|---------------------------|---------------------------------------------------------------|
| removeDeadCellsManualGate | <i>remove dead cells from a flowFrame using manual gating</i> |
|---------------------------|---------------------------------------------------------------|

---

### Description

remove dead cells from a flowFrame, using manual gating in the FSC-A, '(a)Live/Dead' 2D representation. The function uses flowCore::polygonGate()

**Usage**

```
removeDeadCellsManualGate(
  ff,
  preTransform = FALSE,
  transList = NULL,
  FSCChannel,
  LDMarker,
  gateData,
  ...
)
```

**Arguments**

|                           |                                                                                                                                                                                                                                |
|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <code>ff</code>           | a <code>flowCore::flowFrame</code>                                                                                                                                                                                             |
| <code>preTransform</code> | boolean, if TRUE: the <code>transList</code> list of scale transforms will be applied first on the LD channel.                                                                                                                 |
| <code>transList</code>    | applied in conjunction with <code>preTransform == TRUE</code>                                                                                                                                                                  |
| <code>FSCChannel</code>   | a character containing the exact name of the forward scatter channel                                                                                                                                                           |
| <code>LDMarker</code>     | a character containing the exact name of the marker corresponding to (a)Live/Dead channel, or the Live/Dead channel name itself                                                                                                |
| <code>gateData</code>     | a numerical vector containing the polygon gate coordinates first the <code>FSCChannel</code> channel coordinates of each points of the polygon gate, then the LD channel coordinates of each points (prior to scale transform) |
| <code>...</code>          | additional parameters passed to <code>flowCore::polygonGate()</code>                                                                                                                                                           |

**Value**

a `flowCore::flowFrame` with removed dead cells from the input

**Examples**

```
rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
sampleFiles <-
  file.path(rawDataDir, list.files(rawDataDir, pattern = "Donor"))

truncateMaxRange <- FALSE
minLimit <- NULL

# create flowCore::flowSet with all samples of a dataset
fsRaw <- readSampleFiles(
  sampleFiles = sampleFiles,
  whichSamples = "all",
  truncate_max_range = truncateMaxRange,
  min.limit = minLimit)

suppressWarnings(ff_m <- removeMarginsPeacoQC(x = fsRaw[[2]]))
```

```
ff_c <-
  compensateFromMatrix(ff_m,
                      matrixSource = "fcs")

remDeadCellsGateData <- c(0, 0, 250000, 250000,
                        0, 650, 650, 0)

ff_lcells <-
  removeDeadCellsManualGate(ff_c,
                           FSCChannel = "FSC-A",
                           LDMarker = "L/D Aqua - Viability",
                           gateData = remDeadCellsGateData)
```

---

```
removeDebrisManualGate
```

*remove debris from a flowFrame using manual gating*

---

## Description

remove debris from a flowFrame, using manual gating in the FSC-A, SSC-A 2D representation. The function internally uses `flowCore::polygonGate()`

## Usage

```
removeDebrisManualGate(ff, FSCChannel, SSCChannel, gateData, ...)
```

## Arguments

|                         |                                                                                                                                                                                                                           |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <code>ff</code>         | a <code>flowCore::flowFrame</code>                                                                                                                                                                                        |
| <code>FSCChannel</code> | a character containing the exact name of the forward scatter channel                                                                                                                                                      |
| <code>SSCChannel</code> | a character containing the exact name of the side scatter channel                                                                                                                                                         |
| <code>gateData</code>   | a numerical vector containing the polygon gate coordinates first the <code>FSCChannel</code> channel coordinates of each points of the polygon gate, then the <code>SSCChannel</code> channel coordinates of each points. |
| <code>...</code>        | additional parameters passed to <code>flowCore::polygonGate()</code>                                                                                                                                                      |

## Value

a `flowCore::flowFrame` with removed debris events from the input

## Examples

```
rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
sampleFiles <-
  file.path(rawDataDir, list.files(rawDataDir, pattern = "Donor"))
```

```

truncateMaxRange <- FALSE
minLimit <- NULL

# create flowCore::flowSet with all samples of a dataset
fsRaw <- readSampleFiles(
  sampleFiles = sampleFiles,
  whichSamples = "all",
  truncate_max_range = truncateMaxRange,
  min.limit = minLimit)

suppressWarnings(ff_m <- removeMarginsPeacoQC(x = fsRaw[[2]]))

ff_c <-
  compensateFromMatrix(ff_m,
    matrixSource = "fcs")

remDebrisGateData <- c(73615, 110174, 213000, 201000, 126000,
  47679, 260500, 260500, 113000, 35000)

ff_cells <-
  removeDebrisManualGate(ff_c,
    FSCChannel = "FSC-A",
    SSCChannel = "SSC-A",
    gateData = remDebrisGateData)

```

---

```
removeDoubletsCytoPipeline
```

*remove doublets from a flowFrame, using CytoPipeline custom algorithm*

---

## Description

Wrapper around `CytoPipeline::singletGate()`. Can apply the `flowStats` function subsequently on several channel pairs, e.g. (FSC-A, FSC-H) and (SSC-A, SSC-H)

## Usage

```
removeDoubletsCytoPipeline(ff, areaChannels, heightChannels, nmads, ...)
```

## Arguments

|                             |                                                                                       |
|-----------------------------|---------------------------------------------------------------------------------------|
| <code>ff</code>             | a <code>flowCore::flowFrame</code>                                                    |
| <code>areaChannels</code>   | a character vector containing the name of the 'area type' channels one wants to use   |
| <code>heightChannels</code> | a character vector containing the name of the 'height type' channels one wants to use |

`nmads` a numeric vector with the bandwidth above the ratio allowed, per channels pair (cells are kept if the ratio between -A channel[i] and -H channel[i] is smaller than the median ratio + `nmads[i]` times the median absolute deviation of the ratios). Default is 4, for all channel pairs.

`...` additional parameters passed to `CytoPipeline::singletGate()`

### Value

a `flowCore::flowFrame` with removed doublets events from the input

### Examples

```
rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
sampleFiles <-
  file.path(rawDataDir, list.files(rawDataDir, pattern = "Donor"))

truncateMaxRange <- FALSE
minLimit <- NULL

# create flowCore::flowSet with all samples of a dataset
fsRaw <- readSampleFiles(
  sampleFiles = sampleFiles,
  whichSamples = "all",
  truncate_max_range = truncateMaxRange,
  min.limit = minLimit)

suppressWarnings(ff_m <- removeMarginsPeacoQC(x = fsRaw[[2]]))

ff_c <-
  compensateFromMatrix(ff_m,
    matrixSource = "fcs")

ff_s <-
  removeDoubletsCytoPipeline(ff_c,
    areaChannels = c("FSC-A", "SSC-A"),
    heightChannels = c("FSC-H", "SSC-H"),
    nmads = c(3, 5))
```

---

`removeMarginsPeacoQC` *remove margin events using PeacoQC*

---

### Description

Wrapper around `PeacoQC::RemoveMargins()`. Also pre-selects the channels to be handled ( $\Rightarrow$  all signal channels) If input is a `flowSet`, it applies `removeMargins()` to each `flowFrame` of the `flowSet`.

**Usage**

```
removeMarginsPeacoQC(x, channelSpecifications = NULL, ...)
```

**Arguments**

**x** a flowCore::flowSet or a flowCore::flowFrame

**channelSpecifications** A list of lists with parameter specifications for certain channels. This parameter should only be used if the values in the internal parameters description is too strict or wrong for a number or all channels. This should be one list per channel with first a minRange and then a maxRange value. This list should have the channel name found back in colnames(flowCore::exprs(ff)), or the corresponding marker name (found in flowCore::pData(flowCore::description(ff))) . If a channel is not listed in this parameter, its default internal values will be used. The default of this parameter is NULL. If the name of one list is set to AllFluoChannels, then the minRange and maxRange specified there will be taken as default for all fluorescent channels (not scatter)

**...** additional parameters passed to PeacoQC::RemoveMargins()

**Value**

either a flowCore::flowSet or a flowCore::flowFrame depending on the input.

**Examples**

```
rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
sampleFiles <-
  file.path(rawDataDir, list.files(rawDataDir, pattern = "Donor"))

truncateMaxRange <- FALSE
minLimit <- NULL
fsRaw <- readSampleFiles(sampleFiles,
  truncate_max_range = truncateMaxRange,
  min.limit = minLimit)
suppressWarnings(ff_m <- removeMarginsPeacoQC(x = fsRaw[[2]]))
ggplotFilterEvents(ffPre = fsRaw[[2]],
  ffPost = ff_m,
  xChannel = "FSC-A",
  yChannel = "SSC-A")
```

---

resetCellIDs

*reset 'Original\_ID' column in a flowframe*


---

**Description**

: on a flowCore::flowFrame, reset 'Original\_ID' column. This column can be used in plots comparing the events pre and post gating. If the 'Original\_ID' column already exists, the function replaces the existing IDs by the user provided ones. If not, an appendCellID() is called.

**Usage**

```
resetCellIDs(ff, eventIDs = seq_len(nrow(ff)))
```

**Arguments**

`ff` a `flowCore::flowFrame`

`eventIDs` an integer vector containing the values to be set in expression matrix, as Original ID's.

**Value**

new `flowCore::flowFrame` containing the amended (or added) 'Original\_ID' column

**Examples**

```
data(OMIP021Samples)

ff <- appendCellID(OMIP021Samples[[1]])

subsample_ff <- subsample(ff, 100, keepOriginalCellIDs = TRUE)

# re-create a sequence of IDs, ignoring the ones before subsampling
reset_ff <- resetCellIDs(subsample_ff)
```

---

|                 |                                           |
|-----------------|-------------------------------------------|
| runCompensation | <i>compensate with additional options</i> |
|-----------------|-------------------------------------------|

---

**Description**

: this is a simple wrapper around the `flowCore::compensate()` utility, allowing to trigger an update of the fluo channel names with a prefix 'comp-' (as in FlowJo)

**Usage**

```
runCompensation(obj, spillover, updateChannelNames = TRUE)
```

**Arguments**

`obj` a `flowCore::flowFrame` or `flowCore::flowSet`

`spillover` compensation object or spillover matrix or a list of compensation objects

`updateChannelNames` if TRUE, add a 'comp-' prefix to all fluorochrome channels (hence does not impact the columns related to FSC, SSC, or other specific keyword like TIME, Original\_ID, File,...) Default TRUE.

**Value**

a new object with compensated data, and possibly updated column names

**Examples**

```
data(OMIP021Samples)

ff <- OMIP021Samples[[1]]
compMatrix <- flowCore::spillover(ff)$SPILL
ff <- runCompensation(ff,
                      spillover = compMatrix,
                      updateChannelNames = TRUE)
```

---

singletsGate

*Clean doublet events from flow cytometry data*


---

**Description**

will adjust a polygon gate aimed at cleaning doublet events from the flowFrame. The main idea is to use the ratio between the two indicated channel as an indicator and select only the events for which this ratio is 'not too far' from the median ratio. More specifically, the computed ratio is  $ch1/(1+ch2)$ . However, instead of looking at a constant range of this ratio, as is done in `PeacoQC::removeDoublets()`, which leads to a semi-conic gate, we apply a parallelogram shaped gate, by keeping a constant range of channel 2 intensity, based on the target ratio range at the mid value of channel 1.

**Usage**

```
singletsGate(
  ff,
  filterId = "Singlets",
  channel1 = "FSC-A",
  channel2 = "FSC-H",
  nmad = 4,
  verbose = FALSE
)
```

**Arguments**

|                       |                                                                                           |
|-----------------------|-------------------------------------------------------------------------------------------|
| <code>ff</code>       | A <code>flowCore::flowframe</code> that contains flow cytometry data.                     |
| <code>filterId</code> | the name for the filter that is returned                                                  |
| <code>channel1</code> | The first channel that will be used to determine the doublet events. Default is "FSC-A"   |
| <code>channel2</code> | The second channels that will be used to determine the doublet events. Default is "FSC-H" |

|         |                                                                                                                                                                            |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| nmad    | Bandwidth above the ratio allowed (cells are kept if their ratio is smaller than the median ratio + nmad times the median absolute deviation of the ratios). Default is 4. |
| verbose | If set to TRUE, the median ratio and width will be printed. Default is FALSE.                                                                                              |

### Value

This function returns a `flowCore::polygonGate`.

### Examples

```
data(OMIP021Samples)

# simple example with one single singlets gate filter
# FSC-A and FSC-H channels are used by default

mySingletsGate <- singletsGate(OMIP021Samples[[1]], nmad = 3)

selectedSinglets <- flowCore::filter(
  OMIP021Samples[[1]],
  mySingletsGate)

ff_l <- flowCore::Subset(OMIP021Samples[[1]], selectedSinglets)

linRange <- c(0, 250000)

ggplotFilterEvents(
  ffPre = OMIP021Samples[[1]],
  ffPost = ff_l,
  xChannel = "FSC-A", xLinearRange = linRange,
  yChannel = "FSC-H", yLinearRange = linRange)

# application of two singlets gates one after the other

singletsGate1 <- singletsGate(OMIP021Samples[[1]], nmad = 3)
singletsGate2 <- singletsGate(OMIP021Samples[[1]],
  channel1 = "SSC-A",
  channel2 = "SSC-H",
  filterId = "Singlets2")

singletCombinedGate <- singletsGate1 & singletsGate2

selectedSinglets <- flowCore::filter(
  OMIP021Samples[[1]],
  singletCombinedGate)

ff_l <- flowCore::Subset(OMIP021Samples[[1]], selectedSinglets)

ggplotFilterEvents(
  ffPre = OMIP021Samples[[1]],
  ffPost = ff_l,
  xChannel = "FSC-A", xLinearRange = linRange,
```

```

yChannel = "FSC-H", yLinearRange = linRange)

ggplotFilterEvents(
  ffPre = OMIP021Samples[[1]],
  ffPost = ff_1,
  xChannel = "SSC-A", xLinearRange = linRange,
  yChannel = "SSC-H", yLinearRange = linRange)

```

---

|           |                                    |
|-----------|------------------------------------|
| subsample | <i>sub-sampling of a flowFrame</i> |
|-----------|------------------------------------|

---

### Description

: sub-samples a flowFrame with the specified number of samples, without replacement. adds also a column 'Original\_ID' if not already present in flowFrame.

### Usage

```
subsample(ff, nEvents, seed = NULL, keepOriginalCellIDs = TRUE, ...)
```

### Arguments

|                     |                                                                                                                                                                                                                    |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ff                  | a flowCore::flowFrame                                                                                                                                                                                              |
| nEvents             | number of events to be obtained using sub-sampling                                                                                                                                                                 |
| seed                | can be set for reproducibility of event sub-sampling                                                                                                                                                               |
| keepOriginalCellIDs | if TRUE, adds (if not already present) a 'OriginalID' column containing the initial IDs of the cell (from 1 to nrow prior to subsampling). if FALSE, does the same, but takes as IDs (1 to nrow after subsampling) |
| ...                 | additional parameters (currently not used)                                                                                                                                                                         |

### Value

new flowCore::flowFrame with the obtained subset of samples

### Examples

```

data(OMIP021Samples)

# take first sample of dataset, subsample 100 events and create new flowFrame
ff <- subsample(OMIP021Samples[[1]], nEvents = 100)

```

---

|                  |                                                        |
|------------------|--------------------------------------------------------|
| updateMarkerName | <i>update marker name of a given flowFrame channel</i> |
|------------------|--------------------------------------------------------|

---

### Description

: in a flowCore::flowFrame, update the marker name (stored in 'desc' of parameters data) of a given channel. Also update the corresponding keyword in the flowFrame.

### Usage

```
updateMarkerName(ff, channel, newMarkerName)
```

### Arguments

|               |                                                         |
|---------------|---------------------------------------------------------|
| ff            | a flowCore::flowFrame                                   |
| channel       | the channel for which to update the marker name         |
| newMarkerName | the new marker name to be given to the selected channel |

### Value

a new flowCore::flowFrame with the updated marker name

### Examples

```
data(OMIP021Samples)

retFF <- updateMarkerName(OMIP021Samples[[1]],
  channel = "FSC-A",
  newMarkerName = "Fwd Scatter-A")
```

---

|                |                                |
|----------------|--------------------------------|
| writeFlowFrame | <i>write flowFrame to disk</i> |
|----------------|--------------------------------|

---

### Description

wrapper around flowCore::write.FCS() or utils::write.csv that discards any additional parameter passed in (...)

**Usage**

```
writeFlowFrame(
  ff,
  dir = ".",
  useFCSFileName = TRUE,
  prefix = "",
  suffix = "",
  format = c("fcs", "csv"),
  csvUseChannelMarker = TRUE,
  ...
)
```

**Arguments**

|                                  |                                                                                                                                                                                |
|----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <code>ff</code>                  | a <code>flowCore::flowFrame</code>                                                                                                                                             |
| <code>dir</code>                 | an existing directory to store the <code>flowFrame</code> ,                                                                                                                    |
| <code>useFCSFileName</code>      | if <code>TRUE</code> filename used will be based on original fcs filename                                                                                                      |
| <code>prefix</code>              | file name prefix                                                                                                                                                               |
| <code>suffix</code>              | file name suffix                                                                                                                                                               |
| <code>format</code>              | either fcs or csv                                                                                                                                                              |
| <code>csvUseChannelMarker</code> | if <code>TRUE</code> (default), converts the channels to the corresponding marker names (where the Marker is not NA). This setting is only applicable to export in csv format. |
| <code>...</code>                 | other arguments (not used)                                                                                                                                                     |

**Value**

nothing

**Examples**

```
rawDataDir <-
  system.file("extdata", package = "CytoPipeline")
sampleFiles <-
  file.path(rawDataDir, list.files(rawDataDir, pattern = "Donor"))

truncateMaxRange <- FALSE
minLimit <- NULL

# create flowCore::flowSet with all samples of a dataset
res <- readSampleFiles(
  sampleFiles = sampleFiles,
  whichSamples = "all",
  truncate_max_range = truncateMaxRange,
  min.limit = minLimit)

ff_c <- compensateFromMatrix(res[[2]], matrixSource = "fcs")
```

```
outputDir <- base::tempdir()
writeFlowFrame(ff_c,
               dir = outputDir,
               suffix = "_fcs_export",
               format = "csv")
```

# Index

- \* **data**
  - OMIP021Samples, [47](#)
- \* **internal**
  - CytoPipeline, [9](#)
- addProcessingStep
  - (handlingProcessingSteps), [37](#)
- aggregateAndSample, [3](#)
- appendCellID, [4](#)
- applyScaleTransforms, [5](#)
- areFluoCols, [6](#)
- areSignalCols, [6](#)
- as.json.CytoProcessingStep
  - (CytoProcessingStep), [16](#)
- as.list.CytoPipeline
  - (CytoPipeline-class), [10](#)
- as.list.CytoProcessingStep
  - (CytoProcessingStep), [16](#)
- buildCytoPipelineFromCache
  - (interactingWithCytoPipelineCache), [45](#)
- characterOrFunction-class
  - (CytoProcessingStep), [16](#)
- checkCytoPipelineConsistencyWithCache
  - (interactingWithCytoPipelineCache), [45](#)
- cleanProcessingSteps
  - (handlingProcessingSteps), [37](#)
- collectNbOfRetainedEvents
  - (inspectCytoPipelineObjects), [40](#)
- compensateFromMatrix, [7](#)
- computeScatterChannelsLinearScale, [8](#)
- CytoPipeline, [9](#)
- CytoPipeline, character-method
  - (CytoPipeline-class), [10](#)
- CytoPipeline, list-method
  - (CytoPipeline-class), [10](#)
- CytoPipeline, missing-method
  - (CytoPipeline-class), [10](#)
- CytoPipeline-class, [10](#)
- CytoPipeline-class,
  - (CytoPipeline-class), [10](#)
- CytoPipeline-package (CytoPipeline), [9](#)
- CytoPipelineClass, [10](#)
- CytoPipelineClass (CytoPipeline-class), [10](#)
- CytoProcessingStep, [10](#), [16](#)
- CytoProcessingStep-class
  - (CytoProcessingStep), [16](#)
- deleteCytoPipelineCache
  - (interactingWithCytoPipelineCache), [45](#)
- estimateScaleTransforms, [17](#)
- execute, [19](#)
- executeProcessingStep
  - (CytoProcessingStep), [16](#)
- experimentName (CytoPipeline-class), [10](#)
- experimentName<- (CytoPipeline-class), [10](#)
- export2JSONFile (exportCytoPipeline), [24](#)
- exportCytoPipeline, [24](#)
- findTimeChannel, [25](#)
- from.json.CytoProcessingStep
  - (CytoProcessingStep), [16](#)
- getAcquiredCompensationMatrix, [26](#)
- getChannelNamesFromMarkers, [27](#)
- getCPSARGS (CytoProcessingStep), [16](#)
- getCPSFUN (CytoProcessingStep), [16](#)
- getCPSName (CytoProcessingStep), [16](#)
- getCytoPipelineExperimentNames
  - (inspectCytoPipelineObjects), [40](#)

getCytoPipelineFlowFrame  
    (inspectCytoPipelineObjects),  
    40

getCytoPipelineObjectFromCache  
    (inspectCytoPipelineObjects),  
    40

getCytoPipelineObjectInfos  
    (inspectCytoPipelineObjects),  
    40

getCytoPipelineScaleTransform  
    (inspectCytoPipelineObjects),  
    40

getFCSFileName, 28

getNbProcessingSteps  
    (handlingProcessingSteps), 37

getProcessingStep  
    (handlingProcessingSteps), 37

getProcessingStepNames  
    (handlingProcessingSteps), 37

getTransfoParams, 28

ggplotEvents, 30

ggplotFilterEvents, 34

ggplotFlowRate, 36

handlingProcessingSteps, 37

inspectCytoPipelineObjects, 40

interactingWithCytoPipelineCache, 45

OMIP021Samples, 47

pData (CytoPipeline-class), 10

pData<- (CytoPipeline-class), 10

plotCytoPipelineProcessingQueue  
    (inspectCytoPipelineObjects),  
    40

qualityControlFlowAI, 47

qualityControlPeacoQC, 49

readRDSObject, 50

readSampleFiles, 51

removeChannels, 53

removeDeadCellsManualGate, 53

removeDebrisManualGate, 55

removeDoubletsCytoPipeline, 56

removeMarginsPeacoQC, 57

removeProcessingStep  
    (handlingProcessingSteps), 37

resetCellIDs, 58

runCompensation, 59

sampleDisplayNames  
    (CytoPipeline-class), 10

sampleFiles (CytoPipeline-class), 10

sampleFiles<- (CytoPipeline-class), 10

sampleNameFromDisplayName  
    (CytoPipeline-class), 10

show, CytoPipeline-method  
    (CytoPipeline-class), 10

show, CytoProcessingStep-method  
    (CytoProcessingStep), 16

showProcessingSteps  
    (handlingProcessingSteps), 37

singletsGate, 60

subsample, 62

updateMarkerName, 63

writeFlowFrame, 63