## Package 'ddPCRclust'

January 2, 2025

Title Clustering algorithm for ddPCR data

Version 1.26.0

**Description** The ddPCRclust algorithm can automatically quantify the CPDs of nonorthogonal ddPCR reactions with up to four targets. In order to determine the correct droplet count for each target, it is crucial to both identify all clusters and label them correctly based on their position. For more information on what data can be analyzed and how a template needs to be formatted, please check the vignette.

LazyData true

**Depends** R (>= 3.5)

**Imports** plotrix, clue, parallel, ggplot2, openxlsx, R.utils, flowCore, flowDensity (>= 1.13.3), SamSPECTRAL, flowPeaks

Suggests BiocStyle

Collate 'cluster\_functions.R' 'functions.R' 'ddPCRclust.R'

License Artistic-2.0

URL https://github.com/bgbrink/ddPCRclust

BugReports https://github.com/bgbrink/ddPCRclust/issues

biocViews ddPCR, Clustering

RoxygenNote 6.0.1

git\_url https://git.bioconductor.org/packages/ddPCRclust

git\_branch RELEASE\_3\_20

git\_last\_commit 5098ee0

git\_last\_commit\_date 2024-10-29

**Repository** Bioconductor 3.20

Date/Publication 2025-01-02

Author Benedikt G. Brink [aut, cre], Justin Meskas [ctb], Ryan R. Brinkman [ctb]

Maintainer Benedikt G. Brink <bbrink@cebitec.uni-bielefeld.de>

## Contents

calculateCPDs	. 2
createEnsemble	. 3
ddPCRclust	. 4
exportPlots	. 6
exportToCSV	. 7
exportToExcel	. 8
readFiles	
readTemplate	. 9
runDensity	. 10
runPeaks	. 11
runSam	. 12
shearCorrection	. 13
	15

## Index

calculateCPDs (

Calculates the copies per droplet

#### Description

This function takes the results of the clustering and calculates the actual counts per target, as well as the counts per droplet (CPD) for each marker.

#### Usage

```
calculateCPDs(results, template = NULL, constantControl = NULL)
```

#### Arguments

results	The result of the ddPCRclust algorithm.
template	The parsed dataframe containing the template.
constantContro	l
	The constant refrence control, which should be present in each reaction. It is used to normalize the data.

#### Value

A list of lists, containing the counts for empty droplets, each marker with both total droplet count and CPD, and total number of droplets, for each element of the input list respectively.

```
# Read files
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
files <- readFiles(exampleFiles[3])
# To read all example files uncomment the following line
# files <- readFiles(exampleFiles[1:8])
# Read template
```

```
template <- readTemplate(exampleFiles[9])</pre>
```

#### createEnsemble

```
# Run ddPCRclust
result <- ddPCRclust(files, template)</pre>
```

```
# Calculate the CPDs
markerCPDs <- calculateCPDs(result, template$template)</pre>
```

createEnsemble Create a cluster ensemble

#### Description

This function takes the three (or less) clustering approaches of the ddPCRclust package and combines them to one cluster ensemble. See cl\_medoid for more information.

#### Usage

createEnsemble(dens = NULL, sam = NULL, peaks = NULL, file)

#### Arguments

dens	The result of the flowDensity algorithm as a CLUE partition.
sam	The result of the samSPECTRAL algorithm as a CLUE partition.
peaks	The result of the flowPeaks algorithm as a CLUE partition.
file	The input data. More specifically, a data frame with two dimensions, each di- mension representing the intensity for one color.

#### Value

data	The original input data minus the removed events (for plotting)
confidence	The agreement between the different clustering results in percent. If all algo- rithms calculated the same result, the clustering is likely to be correct, thus the confidence is high.
counts	The droplet count for each cluster.

```
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
file <- read.csv(exampleFiles[3])
densResult <- runDensity(file = file, numOfMarkers = 4)
samResult <- runSam(file = file, numOfMarkers = 4)
peaksResult <- runPeaks(file = file, numOfMarkers = 4)
superResult <- createEnsemble(densResult, samResult, peaksResult, file)</pre>
```

ddPCRclust

#### Description

The ddPCRclust algorithm can automatically quantify the events of ddPCR reaction with up to four markers. In order to determine the correct droplet count for each marker, it is crucial to both identify all clusters and label them correctly based on their position. For more information on what data can be analyzed and how a template needs to be formatted, please check the project repository on github.

This is the main function of this package. It automatically runs the ddPCRclust algorithm on one or multiple csv files containing the raw data from a ddPCR run with up to 4 markers.

#### Usage

```
ddPCRclust(files, template, numOfMarkers = 4, sensitivity = 1,
similarityParam = 0.95, distanceParam = 0.2, fast = FALSE,
multithread = FALSE)
```

files	The input data obtained from the csv files. For more information, please see readFiles.	
template	A data frame containing information about the individual ddPCR runs. An example template is provided with this package. For more information, please see readTemplate.	
numOfMarkers	The number of primary clusters that are expected according the experiment set up. Can be ignored if a template is provided. Else, a vector with length equal to length(files) should be provided, containing the number of markers used for the respective reaction.	
sensitivity	A number between 0.1 and 2 determining sensitivity of the initial clustering, e.g. the number of clusters. A higher value means the data is divided into more clusters, a lower value means more clusters are merged. This allows fine tuning of the algorithm for exceptionally low or high CPDs.	
similarityParam		
	If the distance of a droplet between two or more clusters is very similar, it will not be counted for either. The standard it 0.95, i.e. at least 95% similarity. A sensible value lies between 0 and 1, where 0 means none of the 'rain' droplets will be counted and 1 means all droplets will be counted.	
distanceParam	When assigning rain between two clusters, typically the bottom $20\%$ are assigned to the lower cluster and the remaining $80\%$ to the higher cluster. This parameter changes the ratio, i.e. a value of 0.1 would assign only $10\%$ to the lower cluster.	
fast	Run a simpler version of the algorithm that is about 10x faster. For clean data, this might already deliver very good results. However, is is mostly intended to get a quick overview over the data.	
multithread	Distribute the algorithm amongst all CPU cores to speed up the computation.	

#### ddPCRclust

#### Value

results	The results of the ddPCRclust algorithm. It contains three fields:
	data The original input data minus the removed events (for plotting)
	confidence The agreement between the different clustering results in percent
	If all parts of the algorithm calculated the same result, the clustering is likely to
	be correct, thus the confidence is high
	counts The droplet count for each cluster

#### Usage

The main function of the package is ddPCRclust. This function runs the algorithm with one or multiple files, automatically distributing them amongst all cpu cores using the parallel package (parallelization does not work on windows). Afterwards, the results can be exported in different ways, using exportPlots, exportToExcel and exportToCSV. Once the clustering is finished, copies per droplet (CPD) for each marker can be calculated using calculateCPDs.

These functions provide access to all functionalities of the ddPCRclust package. However, expert users can directly call some internal functions of the algorithm, if they find it necessary. Here is a list of all available supplemental functions:

```
runDensity
runSam
runPeaks
createEnsemble
```

#### Author(s)

Maintainer: Benedikt G. Brink <bbrink@cebitec.uni-bielefeld.de>

Other contributors:

- Justin Meskas <jmeskas@bccrc.ca> [contributor]
- Ryan R. Brinkman <rbrinkman@bccrc.ca> [contributor]

#### See Also

Useful links:

- https://github.com/bgbrink/ddPCRclust
- Report bugs at https://github.com/bgbrink/ddPCRclust/issues

```
# Read files
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
files <- readFiles(exampleFiles[3])
# To read all example files uncomment the following line
# files <- readFiles(exampleFiles[1:8])
# Read template
template <- readTemplate(exampleFiles[9])
# Run ddPCRclust
result <- ddPCRclust(files, template)
# Plot the results
```

exportPlots

```
library(ggplot2)
p <- ggplot(data = result$B01$data, mapping = aes(x = Ch2.Amplitude, y = Ch1.Amplitude))
p <- p + geom_point(aes(color = factor(Cluster)), size = .5, na.rm = TRUE) +
ggtitle('B01 example')+theme_bw() + theme(legend.position='none')
p</pre>
```

exportPlots

```
Plot the algorithms results with ggplot2
```

#### Description

A convinience function that takes the results of the ddPCRclust algorithm, plots them using the ggplot2 library and a custom colour palette and saves the plots to a folder.

#### Usage

```
exportPlots(data, directory, annotations, format = "png", invert = FALSE)
```

#### Arguments

data	The result of the ddPCRclust algorithm
directory	The parent directory where the files should saved. A new folder with the exper- iment name will be created (see below).
annotations	Some basic metadata about the ddPCR reaction. If you provided ddPCRclust a template, this paramater can be filled with the corresponding field in the result. Otherwise, you have to provide a character vector containing a name and the the color channels, e.g. c(Name='ddPCR_01-04-2017', Ch1='HEX', Ch2='FAM')
format	Which file format to use. Can be either be a device function (e.g. png), or one of 'eps', 'ps', 'tex' (pictex), 'pdf', 'jpeg', 'tiff', 'png', 'bmp', 'svg' or 'wmf' (windows only). See also ggsave
invert	Invert the axis, e.g. $x = Ch2$ . Amplitude, $y = Ch1$ . Amplitude

#### Value

None

#### Examples

```
# Read files
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
files <- readFiles(exampleFiles[3])
# To read all example files uncomment the following line
# files <- readFiles(exampleFiles[1:8])
# Read template
template <- readTemplate(exampleFiles[9])
# Run ddPCRclust
result <- ddPCRclust(files, template)</pre>
```

# Export the plots

6

#### exportToCSV

```
dir.create('./Results')
exportPlots(data = result, directory = './Results/', annotations = result$annotations)
```

exportToCSV

```
Export the algorithms results to a csv file
```

#### Description

A convinience function that takes the results of the droplClust algorithm and exports them to a csv file.

#### Usage

```
exportToCSV(data, directory, annotations, raw = FALSE)
```

#### Arguments

data	The result of the ddPCRclust algorithm
directory	The parent directory where the files should saved. A new folder with the exper- iment name will be created (see below).
annotations	Some basic metadata about the ddPCR reaction. If you provided ddPCRclust a template, this paramater can be filled with the corresponding field in the result. Otherwise, you have to provide a character vector containing a name and the the color channels, e.g. c(Name='ddPCR_01-04-2017', Ch1='HEX', Ch2='FAM')
raw	Boolean which determines if the annotated raw data should be exported along with the final counts. Basically, a third column will be added to the original data, which contains the cluster number to which this point was assigned to. Useful for example to visualize the clustering later on. (Warning: this can take a while!)

#### Value

None

```
# Read files
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
files <- readFiles(exampleFiles[3])
# To read all example files uncomment the following line
# files <- readFiles(exampleFiles[1:8])
# Read template
template <- readTemplate(exampleFiles[9])
# Run ddPCRclust
result <- ddPCRclust(files, template)
# Export the results
dir.create('./Results')
exportToCSV(data = result, directory = './Results/', annotations = result$annotations)
```

```
exportToExcel
```

#### Description

A convinience function that takes the results of the droplClust algorithm and exports them to an Excel file.

#### Usage

```
exportToExcel(data, directory, annotations, raw = FALSE)
```

#### Arguments

data	The result of the ddPCRclust algorithm
directory	The parent directory where the files should saved. A new folder with the exper- iment name will be created (see below).
annotations	Some basic metadata about the ddPCR reaction. If you provided ddPCRclust a template, this paramater can be filled with the corresponding field in the result. Otherwise, you have to provide a character vector containing a name and the the color channels, e.g. c(Name='ddPCR_01-04-2017', Ch1='HEX', Ch2='FAM')
raw	Boolean which determines if the annotated raw data should be exported along with the final counts. Basically, a third column will be added to the original data, which contains the cluster number to which this point was assigned to. Useful for example to visualize the clustering later on. (Warning: this can take a while!)

#### Value

None

```
# Read files
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
files <- readFiles(exampleFiles[3])
# To read all example files uncomment the following line
# files <- readFiles(exampleFiles[1:8])
# Read template
template <- readTemplate(exampleFiles[9])
# Run ddPCRclust
result <- ddPCRclust(files, template)
# Export the results
dir.create('./Results')
exportToExcel(data = result, directory = './Results/', annotations = result$annotations)
```

readFiles

#### Description

This function reads the raw csv files for ddPCRclust from disk and returns the experiment data. Please refer to the vignette for more information on how these files need to be formatted.

#### Usage

readFiles(files)

#### Arguments

files The input file(s), specifically csv files. Each file represents a two-dimensional data frame. Each row within the data frame represents a single droplet, each column the respective intensities per colour channel.

#### Value

files	A data frame composed of the experiment data
ids	The file ids, e.g. A01, A02, etc.

#### Examples

```
# Read files
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
files <- readFiles(exampleFiles[1:8])</pre>
```

readTemplate

Read a template file from disk

#### Description

This function reads a template file for ddPCRclust from disk and returns a run template and annotations. Please refer to the vignette for information on how this file need to be formatted.

#### Usage

```
readTemplate(template)
```

#### Arguments

template

A csv file containing information about the individual ddPCR runs. An example template is provided with this package. For more information, please check the vignette or the repository on github.

#### Value

annotations	The metatdata provided in the header of the template. It contains four fields:
	Name The name given to this ddPCR experiment
	Ch1 Color channel 1 (usually HEX)
	Ch2 Color channel 2 (usually FAM)
	descriptions Additional descriptions about this ddPCR experiment (e.g. date,
	exprimentor, etc.)
template	A parsed dataframe containing the template.

## Examples

```
# Read template
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
template <- readTemplate(exampleFiles[9])</pre>
```

runDensity

Find the clusters using flowDensity

#### Description

Use the local density function of the flowDensity package to find the cluster centres of the ddPCR reaction. Clusters are then labelled based on their rotated position and lastly the rain is assigned.

#### Usage

```
runDensity(file, sensitivity = 1, numOfMarkers, missingClusters = NULL,
  similarityParam = 0.95, distanceParam = 0.2)
```

file	The input data. More specifically, a data frame with two dimensions, each di- mension representing the intensity for one color channel.	
sensitivity	A number between 0.1 and 2 determining sensitivity of the initial clustering, e.g. the number of clusters. A higher value means more clusters are being found. Standard is 1.	
numOfMarkers	The number of primary clusters that are expected according the experiment set up.	
missingClusters		
	A vector containing the number of primary clusters, which are missing in this dataset according to the template.	
similarityParam		
	If the distance of a droplet between two or more clusters is very similar, it will not be counted for either. The standard it 0.95, i.e. at least 95% similarity. A sensible value lies between 0 and 1, where 0 means none of the 'rain' droplets will be counted and 1 means all droplets will be counted.	
distanceParam	When assigning rain between two clusters, typically the bottom 20% are assigned to the lower cluster and the remaining 80% to the higher cluster. This parameter changes the ratio, i.e. a value of 0.1 would assign only 10% to the lower cluster.	

#### runPeaks

#### Value

data	The original input data minus the removed events (for plotting)
counts	The droplet count for each cluster.
firstClusters	The position of the primary clusters.
partition	The cluster numbers as a CLUE partition (see clue package for more information).

#### Examples

```
# Run the flowDensity based approach
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
file <- read.csv(exampleFiles[3])
densResult <- runDensity(file = file, numOfMarkers = 4)
# Plot the results
library(ggplot2)
p <- ggplot(data = densResult$data, mapping = aes(x = Ch2.Amplitude, y = Ch1.Amplitude))
p <- p + geom_point(aes(color = factor(Cluster)), size = .5, na.rm = TRUE) +
    ggtitle('flowDensity example')+theme_bw() + theme(legend.position='none')
p
```

```
runPeaks
```

Find the clusters using flowPeaks

#### Description

Find the rain and assign it based on the distance to vector lines connecting the cluster centres.

#### Usage

```
runPeaks(file, sensitivity = 1, numOfMarkers, missingClusters = NULL,
  similarityParam = 0.95, distanceParam = 0.2)
```

file	The input data. More specifically, a data frame with two dimensions, each di- mension representing the intensity for one color channel.
sensitivity	A number between 0.1 and 2 determining sensitivity of the initial clustering, e.g. the number of clusters. A higher value means more clusters are being found. Standard is 1.
numOfMarkers	The number of primary clusters that are expected according the experiment set
	up.
missingClusters	
	A vector containing the number of primary clusters, which are missing in this dataset according to the template.
similarityParam	
	If the distance of a droplet between two or more clusters is very similar, it will not be counted for either. The standard it 0.95, i.e. at least 95% similarity. A sensible value lies between 0 and 1, where 0 means none of the 'rain' droplets will be counted and 1 means all droplets will be counted.

distanceParam	When assigning rain between two clusters, typically the bottom 20% are as-
	signed to the lower cluster and the remaining 80% to the higher cluster. This
	parameter changes the ratio, i.e. a value of 0.1 would assign only 10% to the
	lower cluster.

#### Value

data	The original input data minus the removed events (for plotting)
counts	The droplet count for each cluster.
firstClusters	The position of the primary clusters.
partition	The cluster numbers as a CLUE partition (see clue package for more informa- tion).

#### Examples

```
# Run the flowPeaks based approach
exampleFiles <- list.files(paste0(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
file <- read.csv(exampleFiles[3])
peaksResult <- runPeaks(file = file, numOfMarkers = 4)
# Plot the results
library(ggplot2)
p <- ggplot(data = peaksResult$data, mapping = aes(x = Ch2.Amplitude, y = Ch1.Amplitude))
p <- p + geom_point(aes(color = factor(Cluster)), size = .5, na.rm = TRUE) +
ggtitle('flowPeaks example')+theme_bw() + theme(legend.position='none')
p
```

runSam

Find the clusters using SamSPECTRAL

#### Description

Find the rain and assign it based on the distance to vector lines connecting the cluster centres.

## Usage

```
runSam(file, sensitivity = 1, numOfMarkers, missingClusters = NULL,
similarityParam = 0.95, distanceParam = 0.2)
```

file	The input data. More specifically, a data frame with two dimensions, each di- mension representing the intensity for one color channel.
sensitivity	A number between 0.1 and 2 determining sensitivity of the initial clustering, e.g. the number of clusters. A higher value means more clusters are being found. Standard is 1.
numOfMarkers	The number of primary clusters that are expected according the experiment set up.

	missingClusters	
		A vector containing the number of primary clusters, which are missing in this dataset according to the template.
	similarityPara	n
		If the distance of a droplet between two or more clusters is very similar, it will not be counted for either. The standard it 0.95, i.e. at least 95% similarity. A sensible value lies between 0 and 1, where 0 means none of the 'rain' droplets will be counted and 1 means all droplets will be counted.
	distanceParam	When assigning rain between two clusters, typically the bottom $20\%$ are assigned to the lower cluster and the remaining $80\%$ to the higher cluster. This parameter changes the ratio, i.e. a value of 0.1 would assign only $10\%$ to the lower cluster.
Value		
	data	The original input data minus the removed events (for plotting)

uata	The original input data minus the removed events (for plotting)
counts	The droplet count for each cluster.
firstClusters	The position of the primary clusters.
partition	The cluster numbers as a CLUE partition (see clue package for more information).

#### Examples

```
# Run the SamSPECTRAL based approach
exampleFiles <- list.files(paste@(find.package('ddPCRclust'), '/extdata'), full.names = TRUE)
file <- read.csv(exampleFiles[3])
samResult <- runSam(file = file, numOfMarkers = 4)
# Plot the results
library(ggplot2)
p <- ggplot(data = samResult$data, mapping = aes(x = Ch2.Amplitude, y = Ch1.Amplitude))
p <- p + geom_point(aes(color = factor(Cluster)), size = .5, na.rm = TRUE) +
    ggtitle('SamSPECTRAL example')+theme_bw() + theme(legend.position='none')
p
```

shearCorrection Correct for DNA shearing

#### Description

Longer DNA templates produce a lower droplet count due to DNA shearing. This function normalizes the ddPCRclust result based on a stable marker of different lengths to negate the effect of differences in the lengths of the actual markers of interest. (Work in progress)

## Usage

```
shearCorrection(counts, lengthControl, stableControl)
```

## Arguments

counts	The counts per marker as provided by calculateCPDs.
lengthControl	The name of the length Control. If the template name is for example CPT2, the name in the template should be CPT2-125, where 125 represents the number of basepairs.
stableControl	The name of the stable Control used as a reference for this experiment.

## Value

A linear regression model fitting the length vs ln(ratio) (see lm for details on linear regression).

# Index

calculateCPDs, 2, 5, 14 cl\_medoid, 3 createEnsemble, 3, 5 ddPCRclust, 4, 5–8 ddPCRclust-package (ddPCRclust), 4 exportPlots, 5, 6 exportToCSV, 5, 7 exportToExcel, 5, 8 ggsave, 6 lm, <mark>14</mark> parallel, 5 readFiles, 4, 9 readTemplate, 4, 9 runDensity, 5, 10 runPeaks, 5, 11 runSam, 5, 12 shearCorrection, 13